

# BILF150 / BILF200 In-Line Filters

Designed to capture surgical smoke particulate before it enters wall suction systems.



<p><b>SYSTEM 5000®</b> GENERATOR</p>	<p><b>GOLDVAC® SLIM</b> PENCIL</p>	<p><b>GOLDLINE®</b> ESU PENCIL</p>	<p><b>SNAPEVAC™</b> PENCIL ADAPTER</p>

To learn more about these and other innovative products, call **800-448-6506** or visit **CONMED.com**.

## Features & Benefits

Designed to capture surgical smoke particulate before it enters wall suction systems.

- 99.999% ULPA Media
- Economical and easy to use
- Clear housing allows for visualization of filter life



To learn more about this and other innovative products, call **1-800-448-6506** or visit [CONMED.com](http://CONMED.com).

## Ordering Information

	Description	Qty	Catalog Number
<b>NEW</b>	Multi-Use Ideal for connection to wall suction port	10/Case	BILF150
	Multi-Use Ideal for connection to wall canister tree [single 1/4 in (6.4 mm) to 3/8 in (9.5 mm) detached tube]	10/Case	BILF200

CONMED Corporation  
525 French Road, Utica, NY 13502

Customer Service: 1-800-448-6506  
International Sales: 1-315-797-8375  
Fax: 1-800-438-3051

CONMED.com  
info@mail.CONMED.com