




# BioBrace<sup>®</sup>

Reinforced BioInductive Implant



40 mm x 60 mm

## Same Great Benefits, Bigger Size

- Designed for Large Tendon Repair Augmentation: Gluteus Medius, Proximal Hamstring, Extensor Mechanism, Achilles Tendon & More
- Nearly 250% Increase in Surface Area Enables Broad Tissue Coverage
- Provides Up to 650N of Strength<sup>1,3</sup>

## Features & Benefits

- Provides supplemental strength and load-sharing for up to 2 years<sup>1,2</sup>
- Generates new tissue via highly porous collagen scaffold<sup>1,2</sup>
- Fully resorbable<sup>1,2</sup>
- Robust suture compatibility & handling characteristics designed for open or arthroscopic procedures
- Rapidly absorbs 300% of its weight in blood or other fluids
- Off-the-shelf convenience

<sup>1</sup> K203267, K241906, K242187 - 510(k) Clearance Letter - The BioBrace<sup>®</sup> Implant

<sup>2</sup> Based on Preclinical Animal Data

<sup>3</sup> Based on Internal Testing Data

To learn more, call  
**1-866-4CONMED** (Toll Free)  
**727-214-3000** (International)  
or visit [www.CONMED.com](http://www.CONMED.com).



Imagine what we could do together.



## BioBrace® Reinforced BioInductive Implant

Same Great Benefits, Bigger Size

### Ordering Information

Description	Catalog Number
BioBrace® Reinforced BioInductive Implant 40mm x 60mm	BB40X60

### Strengthen Your Repair, **Optimize Healing**



Call 1-866-4CONMED to schedule a trial with your local CONMED Representative. For more information visit us at [www.CONMED.com](http://www.CONMED.com) or scan QR code.

CONMED Corporation  
11311 Concept Blvd.  
Largo, FL 33773

Toll Free: 1-866-4CONMED  
International: 727-214-3000  
[www.CONMED.com](http://www.CONMED.com)

[customerexperience@conmed.com](mailto:customerexperience@conmed.com)  
[internationalorders@conmed.com](mailto:internationalorders@conmed.com)