



# BioBrace<sup>®</sup>

Reinforced BioInductive Implant



5 mm x 250 mm

## Strengthen Your Repair, Optimize Healing

- Designed for Ligament Repair or Reconstruction Augmentation: ACL, Extra Articular Ligaments, UCL, and MPFL
- 42% Increase in Tissue Volume of the ACL When Using BioBrace<sup>®</sup> and a 24% Increase in Maximum Load<sup>1</sup>
- Provides Up To 280N of Supplemental Strength<sup>1,3</sup>

## Features & Benefits

- Provides Supplemental Strength and Load-sharing for Up to 2 years<sup>1,2</sup>
- Generates New Tissue Via Highly Porous Collagen Scaffold<sup>1,2</sup>
- Fully Resorbable<sup>1,2</sup>
- Robust Suture Compatibility & Handling Characteristics Designed for Open or Arthroscopic Procedures
- Rapidly Absorbs 300% of its Weight in Blood or Other Fluids
- Off-the-shelf Convenience

<sup>1</sup> K203267, K241906, K242187 - 510(k) Clearance Letter - The BioBrace<sup>®</sup> Implant

<sup>2</sup> Based on Preclinical Animal Data

<sup>3</sup> Based on Internal Testing Data

To learn more, call  
**1-866-4CONMED** (Toll Free)  
**727-214-3000** (International)  
or visit [www.CONMED.com](http://www.CONMED.com).



Imagine what we could do together.



## BioBrace® Reinforced BioInductive Implant

Strengthen Your Repair, Optimize Healing

### Ordering Information

Description

Catalog Number

BioBrace® Reinforced BioInductive Implant 5mm x 250mm

BB5X250



Call 1-866-4CONMED to schedule a trial with your local CONMED Representative. For more information visit us at [www.CONMED.com](http://www.CONMED.com) or scan QR code.

CONMED Corporation  
11311 Concept Blvd.  
Largo, FL 33773

Toll Free: 1-866-4CONMED  
International: 727-214-3000  
[www.CONMED.com](http://www.CONMED.com)

[customerexperience@conmed.com](mailto:customerexperience@conmed.com)  
[internationalorders@conmed.com](mailto:internationalorders@conmed.com)