

Strengthen Your Repair, Optimize Your Healing

Historically, augment materials have focused on either biologic healing or mechanical strength. BioBrace® is the first and only FDA cleared implant that does both. A highly porous collagen sponge for improved healing reinforced with PLLA microfilaments for strength.

Embrace the Future of Healing

Reinforces

Provides supplemental strength through load sharing ^{1,2}

Generates New Tissue^{1,2}

Biologic scaffold facilitates and optimizes healing ^{1,2}


Resorbs

Supports tissue for ~ 2 years before naturally resorbing ^{1,2}

Why use BioBrace®

- Reinforcement of soft tissue through^{1,2} supplemental strength & load sharing
- Generates new tissue via highly porous 3D DNA absent collagen sponge^{1,2}
- Fully resorbable^{1,2}
- Off-the-shelf convenience
- Suture compatibility and handling characteristics reduce added cost & OR time
- Designed to improve patient outcomes and quicker return to sport/activity

45+ Different Procedures and Counting




¹ K203267 – 510(k) Clearance Letter - The BioBrace® Implant

² Based on preclinical animal data

³ Based on internal test data

Reinforce




BioBrace® has up to
650 newtons³
of Strength

This is equivalent to nearly 150 pounds!


BioBrace Augmented Repair is as Strong as Native Tendon is Just 12 Weeks²

Generates New Tissue



BioBrace® Porous Structure


42% Increase
In ACL Graft Volume²



BioBrace® Filled with New Tissue

180% Increase
In Rotator Cuff Tendon Thickness²

Resorb



BioBrace® protects early repairs and provides supplemental strength for
2 Years
and then fully resorbs¹

Ordering Information

Description	Catalog Number
BioBrace® Reinforced Implant 23mm x 30mm	BB23X30
BioBrace® Reinforced Implant 40mm x 60mm	BB40X60
BioBrace® Reinforced Implant 5mm x 250mm	BB5X250



Strengthen your outcomes today.

Call 1-866-4CONMED to schedule a trial with your local CONMED Representative.
For more information visit us at www.CONMED.com or scan QR code.



Toll Free: 1-866-4CONMED
International: 727-214-3000

customerexperience@conmed.com
internationalorders@conmed.com

www.CONMED.com/BioBrace

©2024 CONMED Corporation, 1363575377 9/24