



## ClearView™ Smoke Evacuation System

CONMED's ClearView™ Smoke Evacuation System is a compact, quiet and effective smoke evacuation solution with state-of-the-art ULPA filtration to combat airborne hazards by removing 99.999% of the smoke contaminants<sup>1</sup>. Surgical smoke contains over 80 chemicals<sup>2,3</sup> that have the same level of mutagenicity as 27-30 cigarettes<sup>4</sup>.



**GOLDVAC® SLIM**  
SMOKE EVACUATION PENCILS



**HELIXAR™ ESU**  
WITH ARGON BEAM COAGULATION



**SYSTEM 5000™**  
ELECTROSURGICAL UNIT



To learn more about these and other innovative products, call 800-448-6506 or visit [CONMED.com](http://CONMED.com)



525 French Road, Utica, NY 13502

CLEARVIEW™ SMOKE  
EVACUATION SYSTEM

## Make your OR smoke-free with CONMED's ClearView™ Smoke Evacuation System.

### Effective Evacuation

System features 4-stage ULPA filtration to remove 99.999% of the smoke contaminants from the surgical site<sup>1</sup>.

### Quiet Operation

Reduces noise pollution in the OR.

### Compact, Easy-to-Integrate Design

Allows for easy placement on shelf, counter or stacking.

### Automatic Activation

System only activates when electro-cautery is used for seamless integration with System 5000 ESU and other ESUs.

### Multiple Settings for Multiple Procedures

Features three flow settings (Low/Lap, Medium & High/Turbo) to meet the needs of different procedure types.

### 3-Port Filter Design

Accommodates a variety of tubing sizes, providing versatility in procedure and accessory setup.

### Up to 35 Hour Filter Life

Delivers unsurpassed value by automatically tracking and recording filter life.

## Ordering Information

Description	Quantity Per Box	Catalog Number
<b>NEW</b> ClearView™ Smoke Evacuator and Accessories		
ClearView™ Smoke Evacuator – 100/120V, Includes: AC Power Cord, Automatic Activation Device, Pneumatic Footswitch, Operator's Manual, and one filter	1 each	60-8585-120
ClearView™ Smoke Evacuator – 220/240V, Includes: AC Power Cord, Automatic Activation Device, Pneumatic Footswitch, Operator's Manual, and one filter	1 each	60-8585-230
ClearView™ Replacement Filter	2/case, 1 pkg	60-8590-001
ClearView™ Automatic Activation Device	1 each	60-8597-001
CONMED Fluid Trap, For use with ClearView™ Smoke Evacuation Unit	10 each	60-8596-001

<sup>1</sup>Author; Environmental Protection Agency EPA-CICA Fact Sheet, EPA-425/F-03-023. Title: Air Pollution Control Technology Fact Sheet [2013]. <sup>2</sup>Al Sahaf OS, Vega-Carrascal I, Cunningham FO, McGrath JP, Bloomfield FJ. Chemical composition of smoke produced by high-frequency electrosurgery. Ir J Med Sci. 2007; 176(3):229-232. <sup>3</sup>Moot AR, Ledingham KM, Wilson PF, et al. Composition of volatile organic compounds in dia thermic plume as detected by selected ion flow tube mass spectrometry. ANZ J Surg. 2007;77 (1-2):20-23. <sup>4</sup>Hill DS, O'Neill JK, Powell RJ, Oliver DW. Surgical smoke—a health hazard in the operating theatre: a study to quantify exposure and a survey of the use of smoke extractor systems in UK plastic surgery units. J Plast Reconstr Aesthet Surg. 2012; 65(7):911-916.

Customer Service: 1-800-448-6506  
International Sales: 1-315-797-8375  
Fax: 1-800-438-3051

CONMED.com  
info@mail.CONMED.com