




DetachaTip[®] II

THE MULTI-USE ENDOSURGERY
INSTRUMENT SOLUTION



An Unmatched Combination of Efficiencies

The DetachaTip® II System from CONMED Corporation provides the ultimate balance between a reusable and disposable hand instrument system. This unique two piece handle and integrated shaft instrument concept delivers excellent cleanability, cost reduction, and consistent performance by enabling you to easily replace high quality, multi-use components as needed. DetachaTip® II eliminates costly repairs and inconsistent performance of conventional reusables. Only DetachaTip® II offers you this combination of efficiencies in performance and affordability.



Performs to the Highest Standards

- Simple two piece, high quality, multi-use components easily assemble and disassemble with no extra parts required.
- Convenient, individual components packaged sterile for immediate introduction into the operative field.
- Integrated shaft/tip assembly delivers consistent, smooth performance for up to 20 sterilization cycles.
- Durable Luer cleaning port permits high-flow flushing of the central shaft for excellent cleaning results.
- Comfortable, durable, ergonomically designed handle provides long-term, repair-free service.
- Eliminates cost, inconvenience, and unreliability of reusable instrument repairs.

SCISSORS



CURVED METZENBAUM



MINI-CURVED METZENBAUM



HOOK

MULTI-USE GRASPERS & DISSECTORS




STRAIGHT FENESTRATED
(DUCKBILL)




ALLIS



5MM BABCOCK




10MM BABCOCK



ENDOWEAVE™



CURVED MARYLAND



TAPERED DOLPHIN

ORDERING & PRODUCT INFORMATION

| DESCRIPTION | QTY | CATALOG NUMBER |
|---|--------|----------------|
| Curved Metzenbaum Scissors, 33cm x 5mm | 1/Unit | 2-1003 |
| Curved Metzenbaum Scissors, 43cm x 5mm | 1/Unit | 2-4301 |
| Curved Mini-Metzenbaum Scissors, 33cm x 5mm | 1/Unit | 2-1004 |
| Hook Scissors, 33cm x 5mm | 1/Unit | 2-1013 |
| Straight Fenestrated (Duckbill) Grasper, 33cm x 5mm | 1/Unit | 2-1008 |
| Straight Fenestrated (Duckbill) Grasper, 43cm x 5mm | 1/Unit | 2-4307 |
| Endoweave™ Grasper, 33cm x 5mm | 1/Unit | 2-1028 |
| Endoweave™ Grasper, 43cm x 5mm | 1/Unit | 2-4328 |
| Allis Grasper, 33cm x 5mm | 1/Unit | 2-1019 |
| Allis Grasper, 43cm x 5mm | 1/Unit | 2-4319 |
| 5mm Babcock Grasper, 33cm x 5mm | 1/Unit | 2-1005 |
| 5mm Babcock Grasper, 43cm x 5mm | 1/Unit | 2-4305 |
| 10mm Babcock Grasper, 33cm x 10mm | 1/Unit | 2-1014 |
| Maryland (Curved) Dissector, 33cm x 5mm | 1/Unit | 2-1009 |
| Maryland (Curved) Dissector, 43cm x 5mm | 1/Unit | 2-4308 |
| Tapered (Dolphin) Dissector, 33cm x 5mm | 1/Unit | 2-1018 |
| Standard Handle with Ratchet | 1/Unit | 1-1010 |
| Luer Port Cap | 5/Box | CD203 |



CONMED Corporation
11311 Concept Blvd.
Largo, FL 33773

Toll Free: 1-866-4CONMED
International: 727-214-3000
www.CONMED.com

customerexperience@conmed.com
internationalorders@conmed.com