

DetachaTip® II

THE MULTI-USE INSTRUMENT SOLUTION

Multi-Use Scissors



CURVED METZENBAUM

Ref. #	Description
2-1003	5mm x 33cm length
2-4301	5mm x 43cm length

- Suggested Uses:**
- General dissection and cutting
 - Spot coagulation



MINI-CURVED METZENBAUM

Ref. #	Description
2-1004	5mm x 33cm length

- Suggested Uses:**
- Nipping the common duct
 - Any precision shearing



HOOK

Ref. #	Description
2-1013	5mm x 33cm length

- Suggested Uses:**
- Nipping the common duct
 - Cutting suture

Multi-Use Graspers & Dissectors



STRAIGHT FENESTRATED (DUCKBILL)

Ref. #	Description
2-1008	5mm x 33cm length
2-4307	5mm x 43cm length

- Suggested Uses:**
- Manipulation of gallbladder, Hartman's pouch, gastric bands
 - GYN applications-very atraumatic



ALLIS

Ref. #	Description
2-1019	5mm x 33cm length
2-4319	5mm x 43cm length

- Suggested Uses:**
- Grasping peritoneal layer of a cyst
 - Manipulation of gallbladder, Hartman's pouch



5MM BABCOCK

Ref. #	Description
2-1005	5mm x 33cm length
2-4305	5mm x 43cm length

- Suggested Uses:**
- Grasping the fundus of the stomach during Nissen Fundoplication
 - Manipulation of bowel



10MM BABCOCK

Ref. #	Description
2-1014	10mm x 33cm length

- Suggested Uses:**
- Grasping the fundus of the stomach during Nissen Fundoplication
 - Manipulation of bowel



ENDOWEAVE™

Ref. #	Description
2-1028	5mm x 33cm length
2-4328	5mm x 43cm length

- Suggested Uses:**
- Counter traction on fundus of gallbladder



CURVED (MARYLAND)

Ref. #	Description
2-1009	5mm x 33cm length
2-4308	5mm x 43cm length

- Suggested Uses:**
- Dissection of common duct and artery
 - Needle passer
 - Hernias—grasping the defect



TAPERED (DOLPHIN)

Ref. #	Description
2-1018	5mm x 33cm length

- Suggested Uses:**
- General dissecting
 - Spot cauterization when closed

Multi-Use Handle



STANDARD

Ref. #	Description
1-1010	Handle with ratchet

