


# ***DURA*glide™**

Stone Removal Balloons



Applied Ingenuity™

**CONMED™**  
ENDOSCOPIC TECHNOLOGIES

# DURAglide™

Stone Removal Balloons

## DURAglide 3™



.035"  
COMPATIBLE


### Triple Lumen Stone Removal Balloons

#### Ordering Information

Catalog No.	Balloon		Catheter		Guidewire Compatibility
	O.D.(mm)	VOL.(cc)	O.D.	Length	
<b>Injection Above Balloon</b>					
007123	8.5	.75	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm
007124	11.5	1.5	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm
007125	15	3.5	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm
<b>Injection Below Balloon</b>					
007143	8.5	.75	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm
007144	11.5	1.5	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm
007145	15	3.5	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm


#### Performance Guaranteed!

##### Triple Lumen



#### Performance Guaranteed!

##### Double Lumen



- Increased Injection Lumen Capacity  
*Increased Flow Reduces Hand Fatigue*
- Color Coded Guidewire Port  
*Visual Indicator Eliminates Confusion*
- 3 Radiopaque 1cm Spaced Distal Marker Bands  
*Enables Fluoroscopic Visualization of Balloon Position*
- Pebax Catheter Shaft  
*Enhances Pushability*



- Standard and Tapered Catheter Options  
*Accommodates a Variety of Duct or Stone Sizes*
- 3 Radiopaque 1cm Spaced Distal Marker Bands  
*Enables Fluoroscopic Visualization of Balloon Position*
- Pebax Catheter Shaft  
*Enhances Pushability*

#### CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road  
Utica, NY 13502  
(315) 797-8375  
Customer Service: 1-800-448-6506  
FAX: 1-800-438-3051  
International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235  
email: info@mail.conmed.com  
www.conmed.com