




# Eliminator<sup>®</sup> PET

Biliary Balloon Dilators



**GORE<sup>®</sup> VIABIL<sup>®</sup>**  
SHORT WIRE BILIARY  
ENDOPROSTHESIS



**FIT<sup>®</sup>**  
BILIARY STONE BALLOON



**TRUPASS<sup>®</sup>**  
SW PAPILLOTOME



To learn more about these and other innovative products, call **1-866-4CONMED** or visit **www.CONMED.com**.




- **Patented High Pressure PET Design**

Does not waist or dogbone

- **Radiopaque Marker Bands**

Provide fluoroscopic identification of balloon working length

- **.035 Guidewire Compatibility**



### Ordering Information

Description	Balloon Size	Catheter OD	Length	Catalog Number
<b>Eliminator® Biliary Balloon Dilators</b>				
Eliminator® PET Biliary Balloon Dilator	12F/4 mm x 2 cm	5F	180 cm	000347
Eliminator® PET Biliary Balloon Dilator	18F/6 mm x 2 cm	5F	180 cm	000348
Eliminator® PET Biliary Balloon Dilator	24F/8 mm x 3 cm	5F	180 cm	000349