



Endoscope Cleaning Brushes

Engineered for optimal cleaning and convenience, CONMED has a broad portfolio of disposable and reusable endoscope cleaning brushes that feature a universal design to meet your clinical needs.



SCOPE SAVER®
BITE BLOCK



CLEARGUARD™
BIOPSY VALVE



OPTIMIZER®
MULTI-CHAMBER POLYP TRAP



To learn more about these and other innovative products, call **1-866-4CONMED** or visit **www.CONMED.com**.

SCRUB SMARTER.

CONMED's Endoscope Cleaning Brushes feature a universal design that's been created to effectively and thoroughly clean a full range of endoscope channels — helping to ensure proper cleaning. All of the brushes in this comprehensive portfolio feature a rounded brush tip to facilitate passage through the scope channel while protecting the biopsy channel from damage.


CHANNEL MASTER™

Disposable Endoscope Cleaning Brushes

Our Channel Master™ series of disposable endoscope cleaning brushes feature a custom crimped shaft designed to effectively traverse tight channel bends and deflect through 90° valve port transitions. To help ensure proper cleaning of the valves, button and biopsy channels, CONMED offers a combination brush featuring a larger head of soft bristles with a convenient "No Slip" handle at one end and channel brush at the other. While the dual ended brush utilizes two brushes with one catheter pass.



Rounded Brush Tip



BLUE BULLET® Disposable Cleaning Brushes

Our complete Blue Bullet® series of disposable cleaning brushes feature blue opaque tubing, a longer working shaft and stiffer bristles to effectively clean all endoscope channels. CONMED offers a combination Blue Bullet® Brush which features a 38mm head with a thumb hold handle to help ensure proper cleaning of the larger valves/channels and 13mm channel brush at the other end.





Endoscopic Cleaning Brushes

Ordering Information

Description	Brush Diameter	Brush Length	Shaft Length	Unit	Catalog Number
-------------	----------------	--------------	--------------	------	----------------

Channel Master™ Endoscope Cleaning Brushes Non-Sterile, Disposable — Biopsy Channel Range 2.0-4.2mm, Shaft O.D. 1.8mm

Pediatric Brush	6.4mm	7.4mm	160cm	50/box	000610
Dual-Ended Brush	6.4mm	7.4mm	230cm	50/box	000615
Single-Ended Brush	6.4mm	7.4mm	230cm	50/box	000639
Combination Valve/Cleaning Brush for Gastroscope/Colonoscope	7.9mm/ 6.4mm	29.2mm/ 7.4mm	230cm	50/box	000617
Large Combination Valve/Cleaning Brush for Gastroscope/Colonoscope	9.7mm/ 6.4mm	45.2mm/ 7.4mm	230cm	50/box	000619*

* Large combination brush includes longer & larger diameter bristles on valve brush. Contains female luer adaptor, hollow shaft, and port hole for fluid injection to soften channel obstructions.

Blue Bullet® Disposable Cleaning Brushes Non-Sterile, Disposable — Biopsy Channel Range 2.0-5.5mm

Bronchoscope	3mm	13mm	115mm	50/box	480
Gastro/Duodeno/Colonoscope	3mm	13mm	250cm	50/box	473
Gastro/Duodeno/Colonoscope	4mm	13mm	250cm	50/box	474
Gastro/Duodeno/Colonoscope	5mm	13mm	250cm	50/box	478
Gastro/Duodeno/Colonoscope	6mm	13mm	250cm	50/box	476
Double End Cleaning Brush for Gastroscope/Colonoscope	5mm	13mm	250cm	50/box	478D
Combination Valve/Cleaning Brush for Gastroscope/Colonoscope	10mm/ 5mm	38mm/ 13mm	250cm	50/box	478C

The 400 series Blue Bullet® feature a blue opaque tubing.

CONMED Corporation
11311 Concept Blvd.
Largo, FL 33773

Toll Free: 1-866-4CONMED
International: 727-214-3000
www.CONMED.com

customerexperience@conmed.com
internationalorders@conmed.com