

EZStart™ Interference Screw

CONMED's NEW EZStart™ Interference Screw is designed to provide strong fixation of bone-tendon-bone (BTB) grafts during ACL reconstruction. This titanium self-tapping, self-notching, interference screw makes initial placement easy – saving steps with an efficient, reproducible technique.



Features and Benefits

Self - Notching

- The tapered nose of the EZStart™ Interference Screw allows for easy initial placement and accurate guidance of the screw between the bone graft and the tunnel.

Self - Tapping

- The innovative dual-start thread design allows EZStart™ Interference Screws to advance when turned, while creating its own thread pattern.
- Designed to eliminate the need for traditional taps.

“ *This self-notching and self-tapping screw saves steps and makes placing in metal screws simple. It bites in the first turn of the driver.*”*

- **Dr. Brett Owens, East Providence, RI**





525 French Road, Utica, NY 13502

EZSTART™
INTERFERENCE
SCREW

Ordering Information

NEW

Description	Dimensions	Catalog Number
EZStart™ Interference Screw		
Interference Screw	6mm x 15mm	EZ0615T
Interference Screw	6mm x 20mm	EZ0620T
Interference Screw	7mm x 15mm - Coming Soon*	EZ0715T
Interference Screw	7mm x 20mm	EZ0720T
Interference Screw	7mm x 25mm	EZ0725T
Interference Screw	7mm x 30mm	EZ0730T
Interference Screw	8mm x 15mm - Coming Soon*	EZ0815T
Interference Screw	8mm x 20mm	EZ0820T
Interference Screw	8mm x 25mm	EZ0825T
Interference Screw	8mm x 30mm	EZ0830T
Interference Screw	9mm x 15mm - Coming Soon*	EZ0915T
Interference Screw	9mm x 20mm	EZ0920T
Interference Screw	9mm x 25mm	EZ0925T
Interference Screw	9mm x 30mm	EZ0930T
Interference Screw	10mm x 15mm - Coming Soon*	EZ1015T
Interference Screw	10mm x 20mm - Coming Soon*	EZ1020T
Interference Screw	10mm x 25mm - Coming Soon*	EZ1025T
Interference Screw	10mm x 30mm - Coming Soon*	EZ1030T
Additional Instrumentation		
Modular Driver Attachment, Straight Hex	3.5mm	D8626
Modular Driver Attachment, Thin Shaft Hex	3.5mm	D8627
Ratcheting Handle	—	D8640
HyperFlex® GuideWire	14" x 0.062" dia.	8572
HyperFlex® GuideWire	9" x 0.062 dia.	8571
Optional Graft Protector	7mm	C8507
Optional Graft Protector	8mm / 9mm	C8509

*510(k) pending not for sale or distribution in the US.

Local 727-392-6464
Toll Free 1-800-237-0169

CONMED.com
customer_service@conmed.com