





NEW

## FIT™ Biliary Stone Balloon

**Designed by endoscopists for endoscopists**, CONMED's New FIT™ (Flexible Intelligent Therapy) SW Stone Balloon adjusts to the duct size to ensure the easy removal of biliary stones. The NEW FIT™ Biliary Stone Balloon inflates to ranges from 9mm up to 22mm reducing the need to change devices.




**TRUPASS™ SW**  
BILIARY PAPILOTOME



**ACE™**  
BILIARY GUIDEWIRE



**GORE® VIABIL®**  
SW BILIARY ENDOPROSTHESIS



To learn more about these and other innovative products, call **800-448-6506** or visit **CONMED.com**.

**NEW**

FIT™ SW  
Stone Balloon



CONMED  
ENDOSCOPIC  
TECHNOLOGIES

## One Balloon-Multiple Sizes

Compatible with short-wire techniques, FIT™ SW Stone Balloons inflate from 9mm – 22mm, so you can now perform effective stone removal with one balloon.

### *One Balloon for Any Size Duct*

- Expandable from 9mm-22mm, the FIT balloon adjusts to the caliber of the bile duct
- Helps provide effective occlusion by adjusting size to the desired duct
- Designed to easily accommodate anatomical variations
- Reduces need to exchange devices mid-procedure

### **Did you know?**

**FIT™ is one of, if not the ONLY, short-wire compatible stone balloons currently available that is NOT made with natural rubber latex.**

# Features and Benefits

## Designed to Reduce Resistance & Bursting

- Rugged yet pliable, adjusts to most duct sizes to match the patient anatomy and easily remove stones
- One of the only stone balloons that allows for better visualization due to its transparent material
- Radiopaque marker band for visibility under fluoroscopy

## Dual-Injection Ports Help Reduce Intraductal Pressure

- Choice of proximal or distal injection helps decrease unnecessary contrast injection into hepatic ducts and can help reduce fluoroscopy time with dual injection ports.

**NEW**

Short-wire  
Compatible



CONMED  
ENDOSCOPIC  
TECHNOLOGIES



INFLATION & INJECTION  
SYRINGES

## *A Complete ERCP Solution:*

*FIT™ SW Stone Balloons complement CONMED's TruPass™ SW Papillotome, ACE™ .025 Biliary Guidewire and GORE® VIABIL® Biliary Endoprosthesis.*



525 French Road, Utica, NY 13502

CONMED  
ENDOSCOPIC  
TECHNOLOGIES

To schedule a demonstration contact your local **CONMED Sales Representative** or, visit **CONMED.com** for more information.

## Ordering Information

Description	Max Vol.	Cath Length	Cath O.D.	Guidewire Compatibility	Catalog Number
<b>NEW</b> FIT™ Biliary Stone Balloon					
FIT™ Short Wire Stone Balloon, 9-22mm	4.5mL	210cm	7F	.035" / .89mm	F922SW
FIT™ Long Wire Stone Balloon, 9-22mm	4.5mL	210cm	7F	.035" / .89mm	F922LW

Not made with natural rubber latex.

Customer Service: 1-800-448-6506  
International Sales: 1-315-797-8375  
Fax: 1-800-438-3051

CONMED.com  
info@mail.CONMED.com