





FlexiTip™

Disposable Sclerotherapy Needle

Flexible Tip, Visible Performance



CLICK-TIP™
INJECTION NEEDLE



BICAP®
BIPOLAR HEMOSTASIS PROBE



CLEARGUARD®
BIOPSY VALVES



To learn more about these and other innovative products, call **800-448-6506** or visit **ConMed.com**.

Optic Yellow Tip

- Optimal visibility in bloody field

Ultra Flexible Tip

- Reliable flexibility for passing over deflected elevator of duodenoscope

Retractable, Lockable Needle

- Minimized needle bounce
- Optimized control during procedure

Ergonomic Handle Design

- Provides control
- Convenient introduction of injection media

OPTIC YELLOW TIP
FOR OPTIMAL VISIBILITY
IN BLOODY FIELD

STANDARD NEEDLE

Ordering Information

Description	Needle Length	Needle Gauge	Sheath Length	Sheath O.D.	Unit	Catalog #
Optic Yellow Tip Needle	4mm	25g	160cm	2.3mm	5/box	000215
Optic Yellow Tip Needle	5mm	25g	160cm	2.3mm	5/box	000216
Optic Yellow Tip Needle	6mm	25g	160cm	2.3mm	5/box	000217
Optic Yellow Tip Colonoscope Needle	5mm	25g	230cm	2.3mm	5/box	000218
Standard Needle	4mm	25g	160cm	2.3mm	5/box	000198
Standard Needle	5mm	25g	160cm	2.3mm	5/box	000199
Standard Needle	6mm	25g	160cm	2.3mm	5/box	000200
Standard Colonoscope Needle	5mm	25g	230cm	2.3mm	5/box	000138

Note: Sterile, Disposable

Customer Service: 1-800-448-6506
International Sales: 1-315-797-8375
Fax: 1-800-438-3051

ConMed.com
info@mail.conmed.com