



# Power Forward with the Hall® Competitive Collection

This series from Hall® includes CONMED saw compatible blades that provide increased cutting efficiency and aggressiveness.



- Premium, Competitive, and Offset Large Bone Oscillating Saw Blades
- Coarse and Fine Tooth Micro Blades
- Large Bone Reciprocating Blades

To learn more, call  
1-866-4CONMED (Toll Free)  
1-727-219-3310 (International)  
or visit [www.CONMED.com](http://www.CONMED.com).

Imagine what we could do together.

# Ordering Information




## Competitive Series

Large bone oscillating saw blades featuring equivalent to better performance compared to Stryker's Dual Cut, Performance, and 2108 series.\*



## Micro Competitive Series

Expanded micro sagittal blade portfolio featuring new coarse and fine-tooth profiles for small bone applications that provide increased cutting efficiency and aggressiveness with a greater election of dimensions.



## Premium Series

The unique honeycomb texture is designed to reduce heat generation and prolonged contact within the cutting block. The tooth profile features teeth with specially ground clearance to help deliver increased cutting efficiency.



## Offset Series

Additional large bone oscillating blades with unique dimensions.



## Performance Cut Series

Designed to be a competitive offering to Stryker's "Dual Cut" blade line.



## Competitive Reciprocator Series

Expanded collection of large bone reciprocating blades.

\*PDD2676205

# Ordering Information

## Choose the blade that fits and convert with confidence.

Use the chart below to find the HALL® compatible blade that aligns with your current dimensions and performance needs. Comparative order codes make it easier to switch from your existing blade while maintaining the size and profile you prefer.

CONMED Saw Compatible Blade Description	CONMED (order codes)	Stryker (comparative order codes)
Competitive Saw Blade Replacement 10 mm X 1.19 mm X 75 mm	CC100-119-75	2108-152-000
Competitive Saw Blade Replacement 10 mm X 1.27 mm X 90 mm	CC100-127-90	4111-127-090
Competitive Saw Blade Replacement 13 mm X 0.89 mm X 90 mm	CC130-089-90	4111-089-090
Competitive Saw Blade Replacement 13 mm X 1 mm X 90 mm	CC130-100-90	2108-103-000
Competitive Saw Blade Replacement 13 mm X 1.19 mm X 90 mm	CC130-119-90	4113-119-090
Competitive Saw Blade Replacement 13 mm X 1.27 mm X 70 mm	CC130-127-70	2108-155-000
Competitive Flared Saw Blade Replacement 13 mm X 1.27 mm X 70 mm	CC130-127-70F	6213-127-070
Competitive Saw Blade Replacement 13 mm X 1.27 mm X 90 mm	CC130-127-90	2108-158-000
Competitive Flared Saw Blade Replacement 17 mm X 1.27 mm X 100 mm	CC170-127-100F	4118-127-100
Competitive Saw Blade Replacement 19.05 mm X 0.89 mm X 90 mm	CC190-089-90	4118-089-075
Competitive Saw Blade Replacement 19.05 mm X 1 mm X 90 mm	CC190-100-90	4118-097-090
Competitive Saw Blade Replacement 19.05 mm X 1.19 mm X 100 mm	CC190-119-100	6118-119-100
Competitive Saw Blade Replacement 19.05 mm X 1.27 mm X 100 mm	CC190-127-100	4118-127-100
Competitive Saw Blade Replacement 19.05 mm X 1.27 mm X 90 mm	CC190-127-90	4118-127-090
Competitive Saw Blade Replacement 22 mm X 1 mm X 90 mm	CC220-100-90	4125-107-090
Competitive Flared Saw Blade Replacement 22 mm X 1.19 mm X 90 mm	CC220-119-90F	6221-119-090
Competitive Saw Blade Replacement 22 mm X 1.27 mm X 90 mm	CC220-127-90	4125-127-090
Competitive Flared Saw Blade Replacement 22 mm X 1.27 mm X 90 mm	CC220-127-90F	6221-127-090
Competitive Flared Saw Blade Replacement 25.4 mm X 0.89 mm X 90 mm	CC254-089-90F	2108-302-000
Competitive Flared Saw Blade Replacement 25.4 mm X 1 mm X 70 mm	CC254-100-70F	4125-089-075
Competitive Flared Saw Blade Replacement 25.4 mm X 1 mm X 90 mm	CC254-100-90F	2108-102-000
Competitive Saw Blade Replacement 25.4 mm X 1.19 mm X 90 mm	CC254-119-90	4125-119-090
Competitive Flared Saw Blade Replacement 25.4 mm X 1.19 mm X 90 mm	CC254-119-90F	6221-119-090
Competitive Saw Blade Replacement 25.4 mm X 1.27 mm X 90 mm	CC254-127-90	4125-127-090
Competitive Flared Saw Blade Replacement 25.4 mm X 1.27 mm X 90 mm	CC254-127-90F	6221-127-090

Imagine what we could do together.

# Ordering Information

CONMED Saw Compatible Blade Description	CONMED (order codes)	Stryker (comparative order codes)
Competitive Flared Saw Blade Replacement 25.4 mm X 1.37 mm X 90 mm	CC254-137-90F	6221-137-090
Competitive Saw Blade Replacement 8 mm X 1.19 mm X 70 mm	CC800-119-70	6213-119-070
Competitive Saw Blade Replacement 8 mm X 1.27 mm X 60 mm	CC800-127-60	2108-127-070
Competitive Precision Micro Blade Replacement 12.5 mm X 0.5 mm X 34.5 mm	CM125-050-35P	2296-003-106
Competitive Micro Blade Replacement 12.5 mm X 0.7 mm X 34.5 mm	CM125-070-35	2296-033-206
Competitive Precision Micro Blade Replacement 13 mm X 0.6 mm X 39 mm	CM130-060-39P	2296-003-108
Competitive Offset Micro Saw Blade Replacement 13 mm X 0.60 mm X 50 mm	CM130-060-50	N/A
Competitive Precision Micro Blade Replacement 17.5 mm X 0.6 mm X 34.5 mm	CM175-060-35P	2296-003-104
Competitive Micro Blade Replacement 5 mm X 0.7 mm X 11.5 mm	CM500-070-12	2296-033-410
Competitive Precision Micro Blade Replacement 5.5 mm X 0.5 mm X 18.5 mm	CM550-050-19P	2296-033-212
Competitive Precision Micro Blade Replacement 5.5 mm X 0.6 mm X 11.5 mm	CM550-060-12P	2296-003-410
Competitive Precision Micro Blade Replacement 5.5 mm X 0.6 mm X 18 mm	CM550-060-18P	2296-003-212
Competitive Precision Micro Blade Replacement 5.5 mm X 0.6 mm X 18.5 mm	CM550-060-19P	2296-003-212
Competitive Precision Micro Blade Replacement 5.5 mm X 0.6 mm X 25 mm	CM550-060-25P	2296-003-414
Competitive Micro Blade Replacement 5.5 mm X 0.7 mm X 18 mm	CM550-070-18	2296-033-212
Competitive Precision Micro Blade Replacement 5.5 mm X 0.7 mm X 18 mm	CM550-070-18P	2296-033-412
Competitive Micro Blade Replacement 5.5 mm X 0.7 mm X 25.5 mm	CM550-070-26	2296-033-414
Competitive Precision Micro Blade Replacement 6.5 mm X 0.6 mm X 18.5 mm	CM650-060-19P	2296-003-051
Competitive Micro Blade Replacement 6.8 mm X 0.6 mm X 18.5 mm	CM680-060-19	2296-003-114
Competitive Precision Micro Blade Replacement 6.8 mm X 0.6 mm X 18.5 mm	CM680-060-19P	2296-033-103
Competitive Micro Blade Replacement 6.8 mm X 0.7 mm X 15 mm	CM680-070-15	2296-033-103
Competitive Precision Micro Blade Replacement 7 mm X 0.6 mm X 29.5 mm	CM700-060-30P	2296-003-115
Competitive Micro Blade Replacement 7 mm X 0.7 mm X 29.5 mm	CM700-070-30	2296-033-115
Competitive Precision Micro Blade Replacement 8 mm X 0.6 mm X 25 mm	CM800-060-25P	2296-003-111
Competitive Micro Blade Replacement 8 mm X 0.6 mm X 30.5 mm	CM800-060-31	2296-003-125
Competitive Micro Blade Replacement 8 mm X 0.7 mm X 25 mm	CM800-070-25	2296-033-111
Competitive Precision Micro Blade Replacement 9.5 mm X 0.6 mm X 11.5 mm	CM950-060-12P	2296-003-101

Imagine what we could do together.

# Ordering Information

CONMED Saw Compatible Blade Description	CONMED (order codes)	Stryker (comparative order codes)
Competitive Precision Micro Blade Replacement 9.5 mm X 0.6 mm X 31 mm	CM950-060-31P	2296-003-225
Competitive Micro Blade Replacement 9.5 mm X 0.7 mm X 31 mm	CM950-070-31	2296-033-225
Competitive Precision Micro Blade Replacement 9.8 mm X 0.6 mm X 18.5 mm	CM980-060-19P	2296-003-105
Competitive Precision Micro Blade Replacement 9.8 mm X 0.6 mm X 25 mm	CM980-060-25P	2296-003-111
Competitive Precision Micro Blade Replacement 9.8 mm X 0.6 mm X 31 mm	CM980-060-31P	2296-003-125
Competitive Micro Blade Replacement 9.8 mm X 0.7 mm X 18.5 mm	CM980-070-19	2296-033-105
Competitive Micro Blade Replacement 9.8 mm X 0.7 mm X 25 mm	CM980-070-25	2296-033-111
Competitive Precision Micro Blade Replacement 9.9 mm X 0.5 mm X 31 mm	CM990-050-31P	2296-003-125
Competitive Offset Saw Blade Replacement 10 mm X 0.8 mm/0.6 mm X 34.5 mm	CO100-080-35	2108-148-000
Competitive Offset Saw Blade Replacement 10.5 mm X 1 mm X 24 mm	CO105-100-24	2108-145-000
Competitive Offset Saw Blade Replacement 10.5 mm X 1 mm X 29.5 mm	CO105-100-29	2108-148-000
Competitive Offset Saw Blade Replacement 11.5 mm X 0.8 mm X 73 mm	CO115-080-73	2108-150-000
Competitive Offset Saw Blade Replacement 12 mm X 1 mm/0.7 mm X 73 mm	CO120-100-73	2108-150-000
Competitive Offset Saw Blade Replacement 18.5 mm X 0.8 mm X 60.5 mm	CO185-080-60	2108-120-000
Competitive Offset Saw Blade Replacement 18.5 mm X 1 mm X 73 mm	CO185-100-73	2108-115-000
Competitive Offset Saw Blade Replacement 22.8 mm X 1 mm X 84 mm	CO228-100-84	2108-105-000
Competitive Offset Saw Blade Replacement 25 mm X 0.8 mm X 73 mm	CO250-080-73	2108-100-000
Competitive Offset Saw Blade Replacement 25 mm X 0.9 mm/0.8 mm X 73 mm	CO250-090-73	2108-105-000
Competitive Offset Saw Blade Replacement 27 mm X 1.2 mm X 62 mm	CO270-120-62	2108-107-001
Competitive Offset Saw Blade Replacement 27.35 mm X 1.4 mm X 66.5 mm	CO2735-140-67	2108-110-000
Competitive Offset Saw Blade Replacement 27.6 mm X 0.8 mm X 58 mm	CO276-080-58	2108-140-000
Competitive Offset Saw Blade Replacement 32 mm X 1 mm X 65 mm	CO320-100-65	2108-107-009
Competitive Offset Saw Blade Replacement 32.5 mm X 1 mm X 57 mm	CO325-100-57	2108-140-000
Competitive Offset Saw Blade Replacement 48 mm X 0.8 mm X 62 mm	CO480-080-62	2108-197-000
Competitive Premium Saw Blade Replacement 13 mm X 1.19 mm X 90 mm	CP130-119-90	6113-119-090
Competitive Premium Saw Blade Replacement 13 mm X 1.27 mm X 90 mm	CP130-127-90	6113-127-090

Imagine what we could do together.

# Ordering Information

CONMED Saw Compatible Blade Description	CONMED (order codes)	Stryker (comparative order codes)
Competitive Premium Saw Blade Replacement 19.05 mm X 1.19 mm X 90 mm	CP190-119-90	6118-119-090
Competitive Premium Saw Blade Replacement 19.05 mm X 1.27 mm X 100 mm	CP190-127-100	6118-127-100
Competitive Premium Saw Blade Replacement 19.05 mm X 1.27 mm X 90 mm	CP190-127-90	6118-127-090
Competitive Premium Saw Blade Replacement 19.05 mm X 1.37 mm X 90 mm	CP190-137-90	6118-137-090
Competitive Premium Saw Blade Replacement 22 mm X 1 mm X 90 mm	CP220-100-90	6121-107-090
Competitive Premium Saw Blade Replacement 22 mm X 1.19 mm X 90 mm	CP220-119-90	6125-119-090
Competitive Premium Flared Saw Blade Replacement 22 mm X 1.19 mm X 90 mm	CP220-119-90F	6221-119-090
Competitive Premium Saw Blade Replacement 22 mm X 1.27 mm X 90 mm	CP220-127-90	6125-127-090
Competitive Premium Flared Saw Blade Replacement 22 mm X 1.27 mm X 90 mm	CP220-127-90F	6221-127-090
Competitive Premium Saw Blade Replacement 22 mm X 1.37 mm X 90 mm	CP220-137-90	6125-137-090
Competitive Premium Saw Blade Replacement 22 mm X 1.47 mm X 90 mm	CP220-147-90	6125-147-090
Competitive Premium Saw Blade Replacement 25.4 mm X 1.27 mm X 100 mm	CP254-127-100	6125-127-100
Competitive Premium Saw Blade Replacement 25.4 mm X 1.27 mm X 90 mm	CP254-127-90	6125-127-090
Competitive Reciprocator Blade Replacement 6.4 mm X 1.14 mm X 60.3 mm	CR064-114-60	0277-096-250
Competitive Reciprocator Blade Replacement 11.7 mm X 1.05 mm X 77.5 mm	CR117-105-78	0277-096-300
Competitive Reciprocator Blade Replacement 11.7 mm X 1.2 mm X 77.5 mm	CR117-120-78	0277-096-326
Competitive Reciprocator Blade Replacement 13.5 mm X 1 mm X 72.5 mm	CR135-100-72	0277-096-275
Competitive Reciprocator Blade Replacement 12.5 mm X 1 mm/.6 mm X 73.5 mm	CR125-106-73	0277-096-275
Competitive Reciprocator Blade Replacement 12.5 mm X 1 mm/.8 mm X 71.1 mm	CR125-108-71	0277-096-276
Competitive Reciprocator Blade Replacement 12.5 mm X 1.2 mm X 73.5 mm	CR125-120-73	0277-096-277
Competitive Precision Micro Premium Blade 5.5 mm X 0.7 mm X 11.5 mm	CM550-070-115P	2296-033-410
Competitive Micro Premium Blade 5.5 mm X 0.6 mm X 18 mm	CM550-060-180	2296-003-212
Competitive Precision Micro Premium Blade 7.0 mm X 0.6 mm X 15 mm	CM700-060-150P	2296-003-103
Competitive Precision Micro Premium Blade 7.0 mm X 0.7 mm X 15 mm	CM700-070-150P	2296-033-103
Competitive Precision Micro Premium Blade 7.0 mm X 0.6 mm X 18.5 mm	CM700-060-185P	2296-003-114
Competitive Precision Micro Premium Blade 7.0 mm X 0.7 mm X 18.5 mm	CM700-070-185P	2296-033-114
Competitive Precision Micro Premium Blade 9.0 mm X 0.6 mm X 11.5 mm	CM900-060-115P	2296-003-411

Imagine what we could do together.


# Ordering Information

CONMED Saw Compatible Blade Description	CONMED (order codes)	Stryker (comparative order codes)
Competitive Precision Micro Premium Blade 9.0 mm X 0.7 mm X 11.5 mm	CM900-070-115P	2296-033-411
Competitive Precision Micro Premium Blade 9.0 mm X 0.6 mm X 18.5 mm	CM900-060-185P	2296-003-105
Competitive Precision Micro Premium Blade 9.0 mm X 0.7 mm X 18.5 mm	CM900-070-185P	2296-033-105
Competitive Precision Micro Premium Blade 9.0 mm X 0.6 mm X 25 mm	CM900-060-250P	2296-003-111
Competitive Precision Micro Premium Blade 9.0 mm X 0.6 mm X 31 mm	CM900-060-310P	2296-003-125
Competitive Precision Micro Premium Blade 9.0 mm X 0.7 mm X 31 mm	CM900-070-310P	2296-033-125
Competitive Precision Micro Premium Blade 13 mm X 0.6 mm X 34.5 mm	CM130-060-345P	2296-003-106
Competitive Precision Micro Premium Blade 13 mm X 0.7 mm X 34.5 mm	CM130-070-345P	2296-003-106
Competitive Precision Micro Premium Blade 13 mm X 0.7 mm X 39 mm	CM130-070-390P	2296-003-108
Competitive Micro Premium Blade 19.5 mm X 0.6 mm X 41 mm	CM195-060-410	2296-003-254
Competitive Precision Micro Premium Blade 19.5 mm X 0.6 mm X 41 mm	CM195-060-413P	2296-003-255
Competitive Micro Premium Blade 17.5 mm X 0.6 mm X 34.5 mm	CM175-060-345	2296-003-504
Competitive Micro Premium Blade 17 mm X 0.7 mm X 34.5 mm	CM170-070-345	2296-033-224
Competitive Micro Premium Blade 20 mm X 0.6 mm X 30 mm	CM200-060-300	5400-134-281
Competitive Reciprocator Blade Replacement 13.5 mm X 1.2 mm X 72.5 mm	CR135-120-72	0277-096-277
Performance Cut Oscillating Saw Blade, 1.27 X 100 X 24,6 mm	DC127-100-25	4125-127-100
Performance Cut Oscillating Saw Blade, 1.27 X 90 X 17,8 mm	DC127-090-18	4118-127-090
Performance Cut Oscillating Saw Blade, 1.19 X 90 X 17,8 mm	DC119-090-18	4118-119-090
Performance Cut Oscillating Saw Blade, 1.19 X 90 X 24,6 mm	DC119-090-25	4125-119-090
Performance Cut Oscillating Saw Blade, 1.27 X 90 X 24,6 mm	DC127-090-25	4125-127-090
Performance Cut Oscillating Saw Blade, 1.27 X 90 X 13 mm	DC127-090-13	4111-127-090
Performance Cut Oscillating Saw Blade, 1.27 X 100 X 17,8 mm	DC127-100-18	4118-127-100
Performance Cut Oscillating Saw Blade, 1.37 X 90 X 13 mm	DC137-090-13	4111-135-090
Performance Cut Oscillating Saw Blade, 1.37 X 90 X 17,8 mm	DC137-090-18	4118-137-090
Performance Cut Oscillating Saw Blade 0.89 mm X 90 mm X 11 mm	DC089-090-11	4111-089-090
Performance Cut Oscillating Saw Blade 0.89 mm X 90 mm X 18 mm	DC089-090-18	4118-089-090
Performance Cut Oscillating Saw Blade 1.19 mm X 85 mm X 21 mm	DC119-085-21	N/A

Imagine what we could do together.

# Ordering Information

CONMED Saw Compatible Blade Description	CONMED (order codes)	Stryker (comparative order codes)
Performance Cut Oscillating Saw Blade 0.89 mm X 75 mm X 25 mm	DC089-075-25	4125-089-075
Performance Cut Oscillating Saw Blade, 0.89 mm X 64 mm X 35 mm	DC089-064-35	4135-089-064
Performance Cut Oscillating Saw Blade 1.27 mm X 64 mm X 35 mm	DC127-064-35	N/A



THIS PRODUCT IS NOT AVAILABLE IN ALL REGIONS

Imagine what we could do together.



CONMED Corporation | 11311 Concept Blvd., Largo, FL 33773 | [www.CONMED.com](http://www.CONMED.com)