



Power Forward with the Hall® Expanded Blade Portfolio.

Introducing the Hall® Competitive Collection

This new series from Hall® includes blades that provide increased cutting efficiency and aggressiveness.

- Premium, Competitive, and Offset Large Bone Oscillating Saw Blades
- Coarse and Fine Tooth Micro Blades
- Large Bone Reciprocating Blades



To learn more, call
1-866-4CONMED (Toll Free)
1-727-219-3310 (International)
or visit www.CONMED.com.

5 NEW BLADE FAMILIES



Competitive Series

Large bone oscillating saw blades featuring equivalent to better performance compared to Stryker's Dual Cut, Performance, and 2108 series.¹



Micro Competitive Series

Expanded micro sagittal blade portfolio featuring new coarse and fine-tooth profiles for small bone applications that provide increased cutting efficiency and aggressiveness with a greater selection of dimensions.



Offset Series

Additional large bone oscillating blades with unique dimensions.



Premium Series

The unique honeycomb texture is designed to reduce heat generation and prolonged contact within the cutting block. The tooth profile features teeth with specially ground clearance to help deliver increased cutting efficiency.



Competitive Reciprocator Series

Expanded collection of large bone reciprocating blades.

Ordering Information



HALL® BLADES

Description	Catalog Number
Competitive Saw Blade Replacement 10mm x 1.19mm x 75mm	CC100-119-75
Competitive Saw Blade Replacement 10mm x 1.27mm x 90mm	CC100-127-90
Competitive Saw Blade Replacement 13mm x 0.89mm x 90mm	CC130-089-90
Competitive Saw Blade Replacement 13mm x 1mm x 90mm	CC130-100-90
Competitive Saw Blade Replacement 13mm x 1.19mm x 90mm	CC130-119-90
Competitive Saw Blade Replacement 13mm x 1.27mm x 70mm	CC130-127-70
Competitive Flared Saw Blade Replacement 13mm x 1.27mm x 70mm	CC130-127-70F
Competitive Saw Blade Replacement 13mm x 1.27mm x 90mm	CC130-127-90
Competitive Flared Saw Blade Replacement 17mm x 1.27mm x 100mm	CC170-127-100F
Competitive Saw Blade Replacement 19.05mm x 0.89mm x 90mm	CC190-089-90
Competitive Saw Blade Replacement 19.05mm x 1mm x 90mm	CC190-100-90
Competitive Saw Blade Replacement 19.05mm x 1.19mm x 100mm	CC190-119-100
Competitive Saw Blade Replacement 19.05mm x 1.27mm x 100mm	CC190-127-100
Competitive Saw Blade Replacement 19.05mm x 1.27mm x 90mm	CC190-127-90
Competitive Saw Blade Replacement 22mm x 1mm x 90mm	CC220-100-90
Competitive Flared Saw Blade Replacement 22mm x 1.19mm x 90mm	CC220-119-90F
Competitive Saw Blade Replacement 22mm x 1.27mm x 90mm	CC220-127-90
Competitive Flared Saw Blade Replacement 22mm x 1.27mm x 90mm	CC220-127-90F
Competitive Flared Saw Blade Replacement 25.4mm x 0.89mm x 90mm	CC254-089-90F
Competitive Flared Saw Blade Replacement 25.4mm x 1mm x 70mm	CC254-100-70F
Competitive Flared Saw Blade Replacement 25.4mm x 1mm x 90mm	CC254-100-90F
Competitive Flared Saw Blade Replacement 25.4mm x 1.19mm x 90mm	CC254-119-90
Competitive Saw Blade Replacement 25.4mm x 1.27mm x 90mm	CC254-119-90F
Competitive Flared Saw Blade Replacement 25.4mm x 1.27mm x 90mm	CC254-127-90
Competitive Flared Saw Blade Replacement 25.4mm x 1.37mm x 90mm	CC254-127-90F
Competitive Saw Blade Replacement 8mm x 1.19mm x 70mm	CC254-137-90F
Competitive Saw Blade Replacement 8mm x 1.27mm x 60mm	CC800-119-70

Ordering Information



HALL® BLADES

Description	Catalog Number
Competitive Saw Blade Replacement 8mm x 1.27mm X 60mm	CC800-127-60
Competitive Precision Micro Blade Replacement 12.5mm x 0.5mm x 34.5mm	CM125-050-35P
Competitive Micro Blade Replacement 12.5mm x 0.7mm x 34.5mm	CM125-070-35
Competitive Precision Micro Blade Replacement 13mm x 0.6mm x 39mm	CM130-060-39P
Competitive Offset Micro Saw Blade Replacement 13mm x 0.60mm X 50mm	CM130-060-50
Competitive Precision Micro Blade Replacement 17.5mm x 0.6mm x 34.5mm	CM175-060-35P
Competitive Micro Blade Replacement 5mm x 0.7mm x 11.5mm	CM500-070-12
Competitive Precision Micro Blade Replacement 5.5mm x 0.5mm x 18.5mm	CM550-050-19P
Competitive Precision Micro Blade Replacement 5.5mm x 0.6mm x 11.5mm	CM550-060-12P
Competitive Precision Micro Blade Replacement 5.5mm x 0.6mm x 18mm	CM550-060-18P
Competitive Precision Micro Blade Replacement 5.5mm x 0.6mm x 18.5mm	CM550-060-19P
Competitive Precision Micro Blade Replacement 5.5mm x 0.6mm x 25mm	CM550-060-25P
Competitive Micro Blade Replacement 5.5mm x 0.7mm x 18mm	CM550-070-18
Competitive Precision Micro Blade Replacement 5.5mm x 0.7mm x 18mm	CM550-070-18P
Competitive Micro Blade Replacement 5.5mm x 0.7mm x 25.5mm	CM550-070-26
Competitive Precision Micro Blade Replacement 6.5mm x 0.6mm x 18.5mm	CM650-060-19P
Competitive Micro Blade Replacement 6.8mm x 0.6mm x 18.5mm	CM680-060-19
Competitive Precision Micro Blade Replacement 6.8mm x 0.6mm x 18.5mm	CM680-060-19P
Competitive Micro Blade Replacement 6.8mm x 0.7mm x 15mm	CM680-070-15
Competitive Precision Micro Blade Replacement 7mm x 0.6mm x 29.5mm	CM700-060-30P
Competitive Micro Blade Replacement 7mm x 0.7mm x 29.5mm	CM700-070-30
Competitive Precision Micro Blade Replacement 8mm x 0.6mm x 25mm	CM800-060-25P
Competitive Micro Blade Replacement 8mm x 0.6mm x 30.5mm	CM800-060-31

Ordering Information



HALL® BLADES

Description	Catalog Number
Competitive Micro Blade Replacement 8mm x 0.7mm x 25mm	CM800-070-25
Competitive Precision Micro Blade Replacement 9.5mm x 0.6mm x 11.5mm	CM950-060-12P
Competitive Precision Micro Blade Replacement 9.5mm x 0.6mm x 31mm	CM950-060-31P
Competitive Micro Blade Replacement 9.5mm x 0.7mm x 31mm	CM950-070-31
Competitive Precision Micro Blade Replacement 9.8mm x 0.6mm x 18.5mm	CM980-060-19P
Competitive Precision Micro Blade Replacement 9.8mm x 0.6mm x 25mm	CM980-060-25P
Competitive Precision Micro Blade Replacement 9.8mm x 0.6mm x 31mm	CM980-060-31P
Competitive Micro Blade Replacement 9.8mm x 0.7mm x 18.5mm	CM980-070-19
Competitive Micro Blade Replacement 9.8mm x 0.7mm x 25mm	CM980-070-25
Competitive Precision Micro Blade Replacement 9.9mm x 0.5mm x 31mm	CM990-050-31P
Competitive Offset Saw Blade Replacement 10mm x 0.8mm/0.6mm x 34.5mm	CO100-080-35
Competitive Offset Saw Blade Replacement 10.5mm x 1mm x 24mm	CO105-100-24
Competitive Offset Saw Blade Replacement 10.5mm x 1mm x 29.5mm	CO105-100-29
Competitive Offset Saw Blade Replacement 11.5mm x 0.8mm x 73mm	CO115-080-73
Competitive Offset Saw Blade Replacement 12mm x 1mm/0.7mm x 73mm	CO120-100-73
Competitive Offset Saw Blade Replacement 18.5mm x 0.8mm x 60.5mm	CO185-080-60
Competitive Offset Saw Blade Replacement 18.5mm x 1mm x 73mm	CO185-100-73
Competitive Offset Saw Blade Replacement 22.8mm x 1mm x 84mm	CO228-100-84
Competitive Offset Saw Blade Replacement 25mm x 0.8mm x 73mm	CO250-080-73
Competitive Offset Saw Blade Replacement 25mm x 0.9mm/0.8mm x 73mm	CO250-090-73
Competitive Offset Saw Blade Replacement 27mm x 1.2mm x 62mm	CO270-120-62
Competitive Offset Saw Blade Replacement 27.35mm x 1.4mm x 66.5mm	CO2735-140-67
Competitive Offset Saw Blade Replacement 27.6mm x 0.8mm x 58mm	CO276-080-58

Ordering Information



HALL® BLADES

Description	Catalog Number
Competitive Offset Saw Blade Replacement 32mm x 1mm x 65mm	CO320-100-65
Competitive Offset Saw Blade Replacement 32.5mm x 1mm x 57mm	CO325-100-57
Competitive Offset Saw Blade Replacement 48mm x 0.8mm x 62mm	CO480-080-62
Competitive Premium Saw Blade Replacement 13mm x 1.19mm x 90mm	CP130-119-90
Competitive Premium Saw Blade Replacement 13mm x 1.27mm x 90mm	CP130-127-90
Competitive Premium Saw Blade Replacement 19.05mm x 1.19mm x 90mm	CP190-119-90
Competitive Premium Saw Blade Replacement 19.05mm x 1.27mm x 100mm	CP190-127-100
Competitive Premium Saw Blade Replacement 19.05mm x 1.27mm x 90mm	CP190-127-90
Competitive Premium Saw Blade Replacement 19.05mm x 1.37mm x 90mm	CP190-137-90
Competitive Premium Saw Blade Replacement 22mm x 1mm x 90mm	CP220-100-90
Competitive Premium Saw Blade Replacement 22mm x 1.19mm x 90mm	CP220-119-90
Competitive Premium Flared Saw Blade Replacement 22mm x 1.19mm x 90mm	CP220-119-90F
Competitive Premium Saw Blade Replacement 22mm x 1.27mm x 90mm	CP220-127-90
Competitive Premium Flared Saw Blade Replacement 22mm x 1.27mm x 90mm	CP220-127-90F
Competitive Premium Saw Blade Replacement 22mm x 1.37mm x 90mm	CP220-137-90
Competitive Premium Saw Blade Replacement 22mm x 1.47mm x 90mm	CP220-147-90
Competitive Premium Saw Blade Replacement 25.4mm x 1.27mm x 100mm	CP254-127-100
Competitive Premium Saw Blade Replacement 25.4mm x 1.27mm x 90mm	CP254-127-90
Competitive Reciprocator Blade Replacement 6.4mm x 1.14mm x 60.3mm	CR064-114-60
Competitive Reciprocator Blade Replacement 11.7mm x 1.05mm x 77.5mm	CR117-105-78
Competitive Reciprocator Blade Replacement 11.7mm x 1.2mm x 77.5mm	CR117-120-78
Competitive Reciprocator Blade Replacement 12.5mm x 1.2mm x 73.5mm	CR125-120-73

Ordering Information



HALL® BLADES

Description	Catalog Number
Competitive Reciprocator Blade Replacement 13.5mm x 1mm x 72.5mm	CR135-100-72
Competitive Reciprocator Blade Replacement 12.5mm x 1mm/.6mm x 73.5mm	CR125-106-73
Competitive Reciprocator Blade Replacement 12.5mm x 1mm/.8mm x 71.1mm	CR125-108-71
Competitive Precision Micro Premium Blade 5.5mm x 0.7mm x 11.5mm	CM550-070-115P
Competitive Precision Micro Premium Blade 5.5mm x 0.6mm x 18mm	CM550-060-180
Competitive Precision Micro Premium Blade 7.0mm x 0.6mm x 15mm	CM700-060-150P
Competitive Precision Micro Premium Blade 7.0mm x 0.7mm x 15mm	CM700-070-150P
Competitive Precision Micro Premium Blade 7.0mm x 0.6mm x 18.5mm	CM700-060-185P
Competitive Precision Micro Premium Blade 7.0mm x 0.7mm x 18.5mm	CM700-070-185P
Competitive Precision Micro Premium Blade 9.0mm x 0.6mm x 11.5mm	CM900-060-115P
Competitive Precision Micro Premium Blade 9.0mm x 0.7mm x 11.5mm	CM900-070-115P
Competitive Precision Micro Premium Blade 9.0mm x 0.6mm x 18.5mm	CM900-060-185P
Competitive Precision Micro Premium Blade 9.0mm x 0.7mm x 18.5mm	CM900-070-185P
Competitive Precision Micro Premium Blade 9.0mm x 0.6mm x 25mm	CM900-060-250P
Competitive Precision Micro Premium Blade 9.0mm x 0.6mm x 31mm	CM900-060-310P
Competitive Precision Micro Premium Blade 9.0mm x 0.7mm x 31mm	CM900-070-310P
Competitive Precision Micro Premium Blade 13mm x 0.6mm x 34.5mm	CM130-060-345P
Competitive Precision Micro Premium Blade 13mm x 0.7mm x 34.5mm	CM130-070-345P
Competitive Precision Micro Premium Blade 13mm x 0.7mm x 39mm	CM130-070-390P
Competitive Micro Premium Blade 19.5mm x 0.6mm x 41mm	CM195-060-410
Competitive Precision Micro Premium Blade 19.5mm x 0.6mm x 41mm	CM195-060-413P
Competitive Micro Premium Blade 17.5mm x 0.6mm x 34.5mm	CM175-060-345



HALL® COMPETITIVE COLLECTION



CONMED Corporation | 11311 Concept Blvd., Largo, FL 33773 | www.CONMED.com



Service & Support