

The XE Factor

The Hall® Titan® family is expanding. The new Titan® XE Modular Drills give you the same reliability and modular compatibility you trust with additional speed options to fit your preferences.

Why Titan XE?

Every surgeon has a preferred feel. The new Titan® XE drills offer optimized reaming speed for those who prefer a slightly slower pace compared to our current Titan® drills, without compromising the power and precision you expect from Hall®.



Key Benefits

- **Expanded Options:** Choose the speed that feels right for you.
- **Same Trusted Design:** Identical look and feel to current Titan® drills now with XE precision.
- **Modular Compatibility:** Works seamlessly with existing Hall® attachments and accessories.
- **Confidence in Performance:** Built on the proven Titan® platform.

Technical Specs at a Glance

Model	Drill Speed	Ream Speed
PRO9100B (Current 1-Trigger)	1500 RPM	300 RPM
PRO9270-1B (New XE 1-Trigger)	1350 RPM	270 RPM
PRO9200B (Current 2-Trigger)	Top Mode: 1500 RPM Bottom Mode: Screw/Tap	Top Mode: 300 RPM Bottom Mode: Screw/Tap
PRO9270-2B (New XE 2-Trigger)	Top Mode: 1500 RPM Bottom Mode: 1350 RPM	Top Mode: 300 RPM Bottom Mode: 270 RPM



Discover the Titan® XE difference.

Contact your sales representative or visit www.CONMED.com to learn more.



Service & Support