



Hi-Fi® Tape

CONMED's widest suture offering, design to maximize tendon compression and prevent tissue tear through.

- 69% less abrasive than FiberTape when measuring tendon tear-through¹
- At 2mm wide, Hi-Fi® Tape offers broader compression than Hi-Fi® Ribbon and #2 Hi-Fi® Suture for increased tendon-to-bone interface
- More than twice as strong as #2 suture²
- Simplifies double-row repairs by eliminating need for medial knot tying

Tendon Tear-Through by Suture Tape



Hi-Fi® Suture

- Braided, non-absorbable suture for soft tissue repair is designed to provide exceptional strength for knot tying while helping minimize tendon damage.
- Peer-reviewed literature has shown that CONMED's #2 Hi-Fi® Suture reduced mean cutting rate by over 50% when compared to FiberWire and FiberTape, which were more abrasive on the suture-tendon interface and led to higher tendon tear-through³.
- Coreless design intended to improve comfort in surgeon's hands when tying knots.
- Available preloaded on all major CONMED soft tissue anchors and individually ranging in size from #2-0 up to #5



Description	Size	Catalog Number
Ribbon		
1.3mm Hi-Fi® Ribbon, 40" Strand, w/Needle (White/Blue)	1.3mm	HR713WN
1.3mm Hi-Fi® Ribbon 40" Strand, no Needle (White/Blue)	1.3mm	HR713W
1.3mm Hi-Fi® Ribbon, 40" Strand, w/Needle (Blue)	1.3mm	HR713BN
1.3mm Hi-Fi® Ribbon, 40" Strand, no Needle (Blue)	1.3mm	HR713B
Also available pre-loaded on Y-Knot® Flex and Y-Knot RC soft tissue anchors		
Tape		
2mm Hi-Fi® Tape (White/Black)	2mm	HT7150
2mm Hi-Fi® Tape (Blue)	2mm	HT7110
Also available pre-loaded Y-Knot® RC and CrossFT® Knotless soft tissue anchors		
Suture		
Hi-Fi® Suture, Single 36" Strand, White-Blue Co-braid HSCX48 0.50in Cutting Needle	#5	H5200
Hi-Fi® Suture, Single 36" Strand, White-Blue Co-Braid, No Needle	#5	H5250
Hi-Fi® Suture, Single 40" Strand, White C-2 0.5in Circle, Tapered Needle	#2	H5000
Hi-Fi® Suture, Single 40" Strand, Blue, 26.00mm 0.50in Circle, Tapered Needle	#2	H6110T
Hi-Fi® Suture, Single 40" Strand, Blue 40.00mm 0.50in Circle, Reverse Cutting Needle	#2	H6140RC
Hi-Fi® Suture, Two 40" Strands, White and White-Blue Co-Braid, C-2 0.50in. Circle, Tapered Needles	#2	H5100
Hi-Fi® Suture, Two 40" Strands, Blue and White-Black Co-Braided, 26.00mm 0.50in Circle Tapered Needles	#2	H6210T
Hi-Fi® Suture, Single 40" Strand, Blue No Needle	#2	H6100
Hi-Fi® Suture, Single 36" Strand, White-Blue Co-Braid, No Needle	#2	H5120
Hi-Fi® Suture, Single 36" Strand, White, No Needle	#2	H5130
Hi-Fi® Suture, Single 36" Strand, White-Black Co-Braid, No Needle	#2	H5150
Hi-Fi® Suture, Single 36" Strand, White-Green Co-Braid, No Needle	#2	H5140
Hi-Fi® Suture, Two 40" Strands Blue and White-Black Co-Braided, No Needle	#2	H6200
Hi-Fi® Suture, Single 36" Strand, White-Blue Co-Braid C-4 0.50in Tapered Needle	#0	H5300
Hi-Fi® Suture, Single 36" Strand, White-Black Co-Braid No Needle	#0	H5301
Hi-Fi® Suture, Single 36" Strand, White, No Needle	#0	H5302
Hi-Fi® Suture, Single 36" Strand, White-Green Co-Braid No Needle	#0	H5303
Hi-Fi® Suture, Single 36" Strand, White-Blue Co-Braid HRT22 0.50in Diamond Point Needle	#0	H5400
Hi-Fi® Suture, Single 36" Strand, White-Blue Co-Braid HR20 0.50in. Tapered Needle	#2-0	H5500
Hi-Fi® Suture, Single 36" Strand, White-Blue Co-Braid HRT20 .50in Diamond Point Needle	#2-0	H5600
Also available pre-loaded Y-Knot RC, Y-Knot Flex, CrossFT®, TruShot™ with Y-Knot soft tissue anchors and more		
Specialty		
Hi-Fi® Passing Loop	#2	HL7000
Hi-Fi® SutureLoop, 60mm Straight Needle, 20" Loop (12 per Box)	#2	HL200
Hi-Fi® SutureLoop, 48mm Curved Needle, 20" Loop (12 per Box)	#2	HL201

CONMED Corporation
11311 Concept Blvd
Largo, FL 33773

Toll Free: 1-866-4CONMED
International: 727-214-3000
www.CONMED.com

customerexperience@conmed.com
internationalorders@conmed.com