



HydroSeal™

Silicone Cannula

Experience seamless access with CONMED's flexible silicone cannula – enhancing maneuverability, safeguarding soft tissues, and maintaining a clear surgical field.

A cannula that manages the workspace so you can focus on the work.

- Low-profile distal flange sits flush against soft tissue and maintains an open workspace
- Large inner diameter facilitates seamless instrument entry & removal
- Double dam manages fluid and minimizes leakage
- Flexible cannula enhances maneuverability inside the workspace
- Removable spacer and ribs along the axis allow for fine-tuning length and the perfect fit in-between sizes



To learn more, call
1-866-4CONMED (Toll Free)
727-214-3000 (International)
or visit www.CONMED.com.

Imagine what we could do together.



HydroSeal™ Silicone Cannula

HydroSeal™ Silicone Cannulas are available in a variety of diameter and length configurations.



Ordering Information

Description	Catalog Number
HydroSeal™ Silicone Cannula 8x20mm	HS-0820
HydroSeal™ Silicone Cannula 8x30mm	HS-0830
HydroSeal™ Silicone Cannula 8x40mm	HS-0840
HydroSeal™ Silicone Cannula 8x50mm	HS-0850
HydroSeal™ Silicone Cannula 10x30mm	HS-1030
HydroSeal™ Silicone Cannula 10x40mm	HS-1040
HydroSeal™ Silicone Cannula 10x50mm	HS-1050
HydroSeal™ Silicone Cannula 12x40mm	HS-1240

Visit www.CONMED.com or connect with a sales representative to learn more.

CONMED Corporation
11311 Concept Blvd.
Largo, FL 33773

Toll Free: 1-866-4CONMED
International: 727-214-3000
www.CONMED.com

customerexperience@conmed.com
internationalorders@conmed.com