



## Specialty Products

From specimen bags and retractors to puncture closure devices and anti-fog kits, CONMED has a comprehensive and expansive portfolio of high-quality, specialty products to meet your clinical and organizational needs.



INDEPENDENT  
MEMORY WIRE BAGS



DISPOSABLE INSUFFLATION  
(VERESS TYPE) NEEDLES



TROCAR WOUND CLOSURE



**AIRSEAL<sup>®</sup> SYSTEM**  
INSUFFLATION MANAGEMENT



**HELIXAR<sup>™</sup> ESU**  
WITH ARGON BEAM COAGULATION



**MICROLAP<sup>®</sup>**  
INSTRUMENTS



To learn more about these and other innovative products, call 800-448-6506 or visit [CONMED.com](http://CONMED.com)

## Disposable Laparoscopic Specimen Bags

Self-opening, high-strength bags with introducer ensure safe specimen capture with easy, one-handed operation. Sterile, single use specimen containers are designed for use in retrieving specimens during endoscopic surgery. Supplied in a dispensing tube for ease of insertion through a 5, 10, or 12mm trocar cannula.

### Unique Design

Patented design provides comfortable ergonomics and can be deployed and manipulated with one hand. Large bags resist tearing, reducing the potential for cell migration.

### Easy to Use

Indicators on the introducer indicate the opening (top) of the bag for correct placement. Bags can be detached within abdominal cavity for later retrieval.

- Minimal insertion force enables smooth deployment
- Insert indicator ensures bag will be in proper position when deployed
- String with slip knot facilitates closure of the specimen bag
- Self-opening bag design provides easy and safe specimen capture
- High strength polyurethane bag material assures no leakage

**NEW**  
Attached Memory  
Wire Bags, Polymer

**NEW**  
Attached Memory  
Wire Bag, Nylon



## Detachable Laparoscopic Specimen Bags

Detachable disposable specimen bags are supplied with an introducer and bag for ease of insertion through a 5, 10, or 12mm trocar cannula.

### Independent Memory Wire Bags

Disposable specimen bags with memory wire in a dispensing tube for ease of insertion through a standard 10, 11, or 12mm cannula.



### Disposable Insufflation (Veress Type) Needles

Used to establish pneumoperitoneum prior to abdominal surgery.

- Retractable safety, spring loaded tip protector. It retracts as the needle is pushed through the abdominal tissue and automatically advances forward once the peritoneum is penetrated.
- Visual needle position indicator
- Clear, comfort handle for easy hang drop test, and a luer lock connector for insufflating the abdominal cavity



### Puncture Closure Devices (Trocar Wound Closure)

Allows full closure of peritoneum and fascial layers through the trocar. Exchangeable through trocar cannula for easy port closure while maintaining pneumoperitoneum. The “J” shaped suture placement needle minimizes risk of organ and tissue injuries. Accommodates both braided and monofilament absorbable sutures, up to size 1.



### Kittner Blunt Dissectors

Convenient, cost-effective packaging houses these sponge and blunt tip probes. Unique dual tip saves time and reduces costs. Non-absorbant radiopaque tip. 40cm length.



### Anti-Fog Kit

Water-based solution reduces fogging of scopes. Includes adhesive backed sponge with radiopaque ribbon.





525 French Road, Utica, NY 13502

## Ordering Information

Description	Qty Per Box	Catalog Number
<b>Anti-Fog Kit</b>		
6cc anti-fog solution vial, sponge single-use, sterile	10/case 1/pkg	CD004
<b>Kittner Blunt Dissectors</b>		
Single tip Kittner, 5mm x 40cm, single use, sterile	5/case 1/pkg	CD801
Dual tip Kittner, 5mm x 40cm, single-use, sterile	5/case 1/pkg	CD802
<b>Luer Port Cap</b>		
Universal Luer Port Cap, non-sterile	5/case 1/pkg	CD203
<b>Puncture Closure Device (Trocar Wound Closure)</b>		
10mm, single-use, sterile, Cannula 10-12mm	10/case 1/pkg	PCD-10
<b>Veress Needle</b>		
Pneumoperitoneum Insufflation Needle 14 gauge, 12cm, with stainless steel tip guard and band indicator	10/case 1/pkg	60-6050-001
Pneumoperitoneum Insufflation Needle 14 gauge, 15cm, with stainless steel tip guard and band indicator	10/case 1/pkg	60-6050-002

Description	Cannula	Diameter	Length	Volume	Qty Per Box	Catalog Number
<b>Specimen Bags</b>						
<b>NEW</b> Attached Memory Wire Bag, Polymer	5mm	3in (7.62cm)	4 in (10.2cm)	120ml	10/case 1/pkg	SB534
Independent Memory Wire Bag, Polymer	10mm	3in (7.62cm)	6in (15.24cm)	255ml	10/case 1/pkg	SB836
Independent Memory Wire Bag, Polymer	10mm	5in (12.7cm)	7in (17.78cm)	680ml	10/case 1/pkg	SB857
Attached Memory Wire Bag, Polymer	10mm	3in (7.62cm)	6in (15.24cm)	255ml	5/case 1/pkg	SB936
Attached Memory Wire Bag, Polymer	10mm	5in (12.7cm)	7in (17.78cm)	680ml	5/case 1/pkg	SB957
Attached Memory Wire Bag, Polymer	12mm	7.5in (19.1cm)	9in (22.9cm)	1500ml	5/case 1/pkg	SB979
<b>NEW</b> Attached Memory Wire Bag, Nylon	10mm	3in (7.62cm)	6in (15.24cm)	255ml	5/case 1/pkg	SB1036
<b>NEW</b> Attached Memory Wire Bag, Nylon	10mm	5in (12.7cm)	7in (17.78cm)	680ml	5/case 1/pkg	SB1057
<b>NEW</b> Attached Memory Wire Bag, Nylon	12mm	7.5in (19.1cm)	9in (22.9cm)	1500ml	5/case 1/pkg	SB1079

Customer Service: 1-800-448-6506  
 International Sales: 1-315-797-8375  
 Fax: 1-800-438-3051

CONMED.com  
 info@mail.CONMED.com