




OptiBite® Disposable Biopsy Forceps


For Your Everyday Biopsy Needs




PRECISOR™ EXL
BIOPSY FORCEP



PRECISOR® HOT
BIOPSY FORCEP



PRECISOR® BRONCHO
BIOPSY FORCEP



To learn more about the innovative CONMED products, call 1-866-4CONMED or visit www.CONMED.com



OVAL CUP



OVAL CUP
W/NEEDLE



ALLIGATOR CUP



ALLIGATOR CUP
W/ NEEDLE



TWIN WIRE
JAW DESIGN

Sheath

- Low friction for optimal pushability
- Reinforced to protect from kinking

Color Coded Cap on Forceps Distinguishes Jaw with and without Needles

- Red cap for jaws with needle
- Yellow cap for jaws without needle
- Protects jaw during handling

Smaller Packaging

- Saves on shelf space



THUMB RING
HANDLE

Ordering Information

Jaw Configuration	Jaw Size	Shaft Length	Unit	Catalog #
Sterile, Coated				
Alligator cup w/needle	2.3mm	230cm	10/box	100575
Alligator cup w/o needle	2.3mm	230cm	10/box	100576
Oval cup w/needle	2.3mm	230cm	10/box	100577
Oval cup w/o needle	2.3mm	230cm	10/box	100578