



CrossFT[®]

Suture Anchors



**Consistent Fixation.
Expanded Options.**

Scan for more information.

Imagine what we could do together.

CrossFT® Anchor

Dual-Thread Fixation for Primary and Medial Row Repair

The CrossFT® platform features PEEK and GENESYS™ biocomposite fully threaded anchors with cortical and cancellous thread profiles designed to provide strong fixation across varying bone qualities.


Fully Threaded Design for Confident Anchor Stability

Built to hold firm and resist migration during repair.



Confident Fixation

The CrossFT® anchor system provides superior fixation strength compared to competitive anchors. ^{1,2,3,4}



CrossFT® GENESYS™ demonstrated **47%** higher pull-out strength compared to competitive systems.^{1,2}
 CrossFT® PEEK demonstrated **20%** higher pull-out strength compared to competitive systems.^{1,2,3}

Imagine what we could do together.

Options for Every Preference


A single anchor platform offering PEEK and GENESYS™ biocomposite materials with double or triple-loaded #2 Hi-Fi® suture and 1.3 mm Hi-Fi® Ribbon options, enabling tailored constructs without changing technique.



CrossFT® with needles for open or mini-open procedures



GENESYS™ CrossFT®
Biocomposite material with #2 Hi-Fi® suture and 1.3 mm Hi-Fi® Ribbon



GENESYS™ CrossFT®
anchors maintain mechanical strength during the critical early healing period.⁵



CrossFT® PEEK radiolucent, non-resorbable polymer with #2 Hi-Fi® suture

Imagine what we could do together.

Consistent Fixation. Expanded Options.

GENESYS™ CrossFT® with Hi-Fi®
Ribbon for broader compression and
expanded procedural flexibility.

Broader Compression

At twice the width of #2 suture, 1.3 mm Hi-Fi® Ribbon is designed to increase tendon-to-bone contact, improve load distribution, and support smooth sliding knot performance.

Low-Profile Knots

Hi-Fi® Ribbon creates smaller knot stacks than #2 suture, helping reduce abrasion and post-operative discomfort associated with larger knots.⁶

Strength You Can Trust

Hi-Fi® Ribbon is 14% stronger than #2 suture for added confidence and control, while GENESYS™ material with microTCP™ is engineered to enhance anchor strength by up to 51%.^{7,8}

Imagine what we could do together.



CrossFT® Tap



CrossFT® Punch



CrossFT®
Universal Punch⁸

Streamlined Instrumentation

Broaches and taps available in
4.5, 5.5, and 6.5 mm sizes with
color-coded handles for
easy identification.



Imagine what we could do together.

Ordering Information

GENESYS™ CrossFT® Suture Anchors	SKU
4.5 mm GENESYS™ CrossFT® Suture Anchor w/ two 1.3 mm Hi-Fi® Ribbons	CFBC-4502R
4.5 mm GENESYS™ CrossFT® Suture Anchor w/ three 1.3 mm Hi-Fi® Ribbons	CFBC-4503R
5.5 mm GENESYS™ CrossFT® Suture Anchor w/ two 1.3 mm Hi-Fi® Ribbons	CFBC-5502R
5.5 mm GENESYS™ CrossFT® Suture Anchor w/ three 1.3 mm Hi-Fi® Ribbons	CFBC-5503R
6.5 mm GENESYS™ CrossFT® Suture Anchor w/ two 1.3 mm Hi-Fi® Ribbons	CFBC-6502R
6.5 mm GENESYS™ CrossFT® Suture Anchor w/ three 1.3 mm Hi-Fi® Ribbons	CFBC-6503R
4.5 mm GENESYS™ CrossFT® w/ two #2 Hi-Fi® Sutures (White, WhiteBlue)	CFBC-4502
4.5 mm GENESYS™ CrossFT® w/ two #2 Hi-Fi® Sutures (Blue, WhiteBlue)	CFBC-4502B
4.5 mm GENESYS™ CrossFT® w/ three #2 Hi-Fi® Sutures (White, WhiteBlue, WhiteBlack)	CFBC-4503
4.5 mm GENESYS™ CrossFT® w/ three #2 Hi-Fi® Sutures (Blue, WhiteBlue, WhiteBlack)	CFBC-4503B
5.5 mm GENESYS™ CrossFT® w/ two #2 Hi-Fi® Sutures (White, WhiteBlue)	CFBC-5502
5.5 mm GENESYS™ CrossFT® w/ two #2 Hi-Fi® Sutures (Blue, WhiteBlue)	CFBC-5502B
5.5 mm GENESYS™ CrossFT® w/ three #2 Hi-Fi® Sutures (White, WhiteBlue, WhiteBlack)	CFBC-5503
5.5 mm GENESYS™ CrossFT® w/ three #2 Hi-Fi® Sutures (Blue, WhiteBlue, WhiteBlack)	CFBC-5503B
6.5 mm GENESYS™ CrossFT® w/ two #2 Hi-Fi® Sutures (White, WhiteBlue)	CFBC-6502
6.5 mm GENESYS™ CrossFT® w/ two #2 Hi-Fi® Sutures (Blue, WhiteBlue)	CFBC-6502B
6.5 mm GENESYS™ CrossFT® w/ three #2 Hi-Fi® Sutures (White, WhiteBlue, WhiteBlack)	CFBC-6503
6.5 mm GENESYS™ CrossFT® w/ three #2 Hi-Fi® Sutures (Blue, WhiteBlue, WhiteBlack)	CFBC-6503B



Imagine what we could do together.

Ordering Information

GENESYS™ CrossFT® Suture Anchors w/ Needles	SKU
4.5 mm GENESYS™ CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFBC-4502N
4.5 mm GENESYS™ CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (Blue, WhiteBlack) (17 mm Length) w/ needles	CFBC-4502NB
4.5 mm GENESYS™ CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (White, WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFBC-4503N
4.5 mm GENESYS™ CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (Blue, WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFBC-4503NB
5.5 mm GENESYS™ CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFBC-5502N
5.5 mm GENESYS™ CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (Blue, WhiteBlack) (17 mm Length) w/ needles	CFBC-5502NB
5.5 mm GENESYS™ CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (White, WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFBC-5503N
5.5 mm GENESYS™ CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (Blue, WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFBC-5503NB
6.5 mm GENESYS™ CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFBC-6502N
6.5 mm GENESYS™ CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (Blue, WhiteBlack) (17 mm Length) w/ needles	CFBC-6502NB
6.5 mm GENESYS™ CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (White, WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFBC-6503N
6.5 mm GENESYS™ CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (Blue, WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFBC-6503NB
GENESYS™ CrossFT® Instruments	SKU
4.5 mm GENESYS™ CrossFT® Broaching Punch	PFT-45MB
5.5 mm GENESYS™ CrossFT® Broaching Punch	PFT-55MB
6.5 mm GENESYS™ CrossFT® Broaching Punch	PFT-65MB
4.5 mm GENESYS™ CrossFT® Tap	TFT-45MB
5.5 mm GENESYS™ CrossFT® Tap	TFT-55MB
6.5 mm GENESYS™ CrossFT® Tap	TFT-65MB



Imagine what we could do together.


Ordering Information

PEEK CrossFT® Suture Anchors	SKU
4.5 mm CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (White, WhiteBlue)	CFP-4502
4.5 mm CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (Blue, WhiteBlue)	CFP-4502B
4.5 mm CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (White, WhiteBlue, WhiteBlack)	CFP-4503
4.5 mm CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (Blue, WhiteBlue, WhiteBlack)	CFP-4503B
5.5 mm CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (White, WhiteBlue)	CFP-5502
5.5 mm CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (Blue, WhiteBlue)	CFP-5502B
5.5 mm CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (White, WhiteBlue, WhiteBlack)	CFP-5503
5.5 mm CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (Blue, WhiteBlue, WhiteBlack)	CFP-5503B
6.5 mm CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (White, WhiteBlue)	CFP-6502
6.5 mm CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (Blue, WhiteBlue)	CFP-6502B
6.5 mm CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (White, WhiteBlue, WhiteBlack)	CFP-6503
6.5 mm CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (Blue, WhiteBlue, WhiteBlack)	CFP-6503B
PEEK CrossFT® Suture Anchors w/ Needles	SKU
4.5 mm CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFP-4502N
4.5 mm CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (Blue, WhiteBlack) (17 mm Length) w/ needles	CFP-4502NB
4.5 mm CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (White, WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFP-4503N
4.5 mm CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (Blue, WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFP-4503NB
5.5 mm CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFP-5502N
5.5 mm CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (Blue, WhiteBlack) (17 mm Length) w/ needles	CFP-5502NB
5.5 mm CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (White, WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFP-5503N
5.5 mm CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (Blue, WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFP-5503NB
6.5 mm CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFP-6502N
6.5 mm CrossFT® Suture Anchor w/ two #2 Hi-Fi® Sutures (Blue, WhiteBlack) (17 mm Length) w/ needles	CFP-6502NB
6.5 mm CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (White, WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFP-6503N
6.5 mm CrossFT® Suture Anchor w/ three #2 Hi-Fi® Sutures (Blue, WhiteBlack, WhiteGreen) (17 mm Length) w/ needles	CFP-6503NB

Imagine what we could do together.

Ordering Information

PEEK CrossFT® Instruments	SKU
Universal CrossFT® Broaching Punch	PFT-00M
4.5 mm CrossFT® Broaching Punch	PFT-45M
6.5 mm GENESYS™ CrossFT® Broaching Punch	PFT-55M
6.5 mm CrossFT® Broaching Punch	PFT-65M
4.5 mm CrossFT® Tap	TFT-45M
5.5 mm CrossFT® Tap	TFT-55M
6.5 mm CrossFT® Tap	TFT-65M



1. Barber FA, et al. Cyclic Loading Biomechanical Analysis of the Pullout Strengths of Rotator Cuff and Glenoid Anchors: 2013 Update. Arthroscopy 2013; 29:832–844.
2. Tashjian RZ, et al. Initial Fixation Strength of Massive Rotator Cuff Tears: In Vitro Comparison of Single-Row Suture Anchor and Transosseous Tunnel Constructs. Arthroscopy 2007; 23:710–716.
3. Barber FA, et al. Biomechanical Analysis of Pullout Strengths of Rotator Cuff and Glenoid Anchors: 2011 Update. Arthroscopy 2011; 27:895-905.
4. TR09-378
5. Data on File; PDD2940135
6. Data on File; TR17-01219
7. Data on File; TR18-00183-1, TR12-441
8. Patent-pending broaching punch is designed to minimize stress fractures.

Imagine what we could do together.



Imagine what we could do together.

CONMED Corporation | 11311 Concept Blvd., Largo, FL 33773 | www.CONMED.com



Service & Support