



OrthoVAC™

Single-Use Large Bore Suction Kit

For use in knee and hip arthroplasty

Big procedures require big tools. OrthoVAC™ offers one of the largest bore Yankauers on the market for reliable, efficient particulate removal.



More Suction, Less Occlusion

OrthoVAC™ is designed to reduce clogging and improve efficiency in your cases. The 6mm Yankauer and 1 1/32" tubing remove large bone fragments and tissue particulate effectively.

All-In-One Customizable Kit

Meeting procedural needs and personal preferences, this single-SKU kit allows surgeons a customizable experience. Kit includes:

- Large Bore Yankauer (1)
 - 1 1/32" tubing
 - Suction canister adapter, blue (1)
 - Standard instrument adapter, green (1)
 - Pour Spout Adapter (1)
- *Kit does not include Canister*

Adapter Versatility

"Blue" and "Pour Spout" adapters provide seamless connectivity, regardless of your facility's suction canister or vacuum system. The "Green" adapter allows suction instruments to be interchangeable.

Committed to Eliminating Cross Contamination

All OrthoVAC™ suction instruments are sterile, single-use devices, designed to increase patient safety and eliminate risk of HAI.

▼ **Don't let a clog press pause on your procedure.**
Visit www.CONMED.com/OrthoVAC to schedule an evaluation today.

Introducing Next Generation Innovation

If you're doing total joints, you need solutions designed to support better outcomes and busy caseloads. That's where we come in.

Hall® Titan® Power Series

The next generation of powerful, precise, and reliable orthopedic power tools.



Hall® Titan®
Modular Drill



Hall® Titan®
Primecut+™ Oscillating Saw



AssistArm® Limb Positioner

Unparalleled freedom when maneuvering - reliable stability when you're not.



Ordering Information

Description	Quantity	Catalog Number
OrthoVAC™ Orthopedic Yankauer Kit (Sterile)	10	0037001

CONMED Corporation
11311 Concept Blvd.,
Largo, FL 33773

Toll Free: 1-866-4CONMED
International: 727-214-3000
www.CONMED.com

customerexperience@CONMED.com
internationalorders@CONMED.com