



PADPRO[®]

Adult/Child & Infant Multifunction
Electrodes for Defibrillation/Pacing/
Cardioversion/Monitoring



- Constructed to adhere and conform to patient's body, providing excellent current distribution
- Adult/child electrodes offer a maximum of 50 shocks, 8 hours of continuous pacing, or 24 hours of patient monitoring
- Pre-connect packaging options available
- 36-month shelf life
- PadPro[®] adapters allow multiple brands of defibrillators to use the same electrodes



ADULT
ECG ELECTRODES



PEDIATRIC
ECG ELECTRODES



PULSE OXIMETRY
SENSORS



To learn more about these and
other innovative products,
call **1.866.4CONMED**
or visit **CONMED.com**.



Ordering Information

Adult/Child & Infant Multifunction Electrodes for Defibrillation/Pacing/Cardioversion/Monitoring



Size A



Size A/B



Size C



Size D



Pre-Connect Packaging

Product Description	Background Label Color	Size	Leadwire Length	Connector Style and Product Label Color		
				Universal/ Anderson ¹	Zoll [®]	Physio-Control [®] Quik-Combo [™]
Radiotransparent, Adult/Child	White	A	60"	2001	2001Z	2001M
Radiotransparent, Adult/Child w/carbon leadwires	White	A	60"	2001-C	2001Z-C	2001M-C
Radiotransparent, Adult/Child w/Pre-Connect Packaging	White	A	60"	2001-PC	2001Z-PC	2001M-PC
Radiotransparent, Adult w/large-posterior pad	White	A/B	60"	2001-EPS	N/A	NA
Radiotransparent, Adult/Child w/Sterile Packaging	White	A	60"	2502	2502Z	2502M
Radiotranslucent, Adult/Child	White	A	60"	2516	2516Z	2516M
Radiotranslucent, Adult/Child w/Pre-Connect Packaging	White	A	60"	2516-PC	2516Z-PC	2516M-PC
Radiotranslucent, Infant	Pink	D	60"	2603	2603Z	2603M
Radiotranslucent, Mini-Infant	Blue	C	24"	2602	2602Z	2602M

Sizes:

A = 6.1"	(15.5cm) x 4.5"	(11.4cm)
B = 9.0"	(22.9cm) x 6.0"	(15.2cm)
C = 2.6"	(6.7cm) x 2.2"	(5.6cm)
D = 4.3"	(11.0cm) x 2.9"	(7.3cm)

¹Electrodes with the Universal/Anderson connector can be utilized in environments with multiple brands of defibrillators (via PadPro[®] adapters).

Product and company names listed are trademarks or trade names of their respective companies.