



PlumePen® Elite Surgical Smoke Evacuation Pencil

From our Buffalo Filter® line of Smoke Evacuation products, PlumePen® Elite features a compact, ultra slim, ergonomic design that is significantly smaller than competitive electrosurgical smoke pencils.



SYSTEM 5000
GENERATOR



GOLDVAC® SLIM
PENCILS



VIROSAFE® VISICLEAR®
FILTER



To learn more about these and other innovative products, call 800-448-6506 or visit CONMED.com.

PlumePen® Elite

Designed for Surgeons by Surgeons

Unrivaled Profile

- Powerfully compact and ultra slim
- Similar in size to standard electro-surgical pencils
- Significantly smaller than competitive electro-surgical smoke pencils



Over-molded Button Design

One piece over-molded button design limits fluid ingress aiding in the prevention of sticking buttons

Full Blade Reveal

Provides exceptional visibility at surgical site

Adjustable Capture Port

Offers customized placement for effective smoke plume capture regardless of blade length

360 Degree Swivel

Allows for fluid hand movement and free range of motion

Unsurpassed Comfort & Performance

Over-molded Cradle Design

- Provides improved grip and comfort

Up to 25% More Flow

- Provides more flow than competitive electro-surgical smoke pencils

Light, Flexible Tubing

- Aids in reducing pull, hand fatigue, and tension on the wrist

Wire Management

- Seamless integration of ESU wire and tubing
- Provides organized cord management on the sterile field

Simple Connectivity

- Easily connects to most surgical plume evacuators and ESU generators



PLUMEPEN® ELITE
SURGICAL SMOKE
EVACUATION PENCIL

Ordering Information

Description	Qty	Catalog Number
PlumePen® Elite (Stainless Steel ESU Blade)		
10 ft. (3 m) tubing	5/Case	PLP1005
10 ft. (3 m) tubing	20/Case	PLP1020
15 ft. (4.6 m) tubing	5/Case	PLP1505
15 ft. (4.6 m) tubing	20/Case	PLP1520
PlumePen® Elite (Non-Stick ESU Blade)		
10 ft. (3 m) tubing	5/Case	PLP2005
10 ft. (3 m) tubing	20/Case	PLP2020
15 ft. (4.6 m) tubing	5/Case	PLP2505
15 ft. (4.6 m) tubing	20/Case	PLP2520

CONMED Corporation
525 French Road, Utica, NY 13502

Local 727-392-6464
Toll Free 1-800-237-0169

CONMED.com
customer_service@CONMED.com