

PlumePort® ActiV®

Laparoscopic Smoke Filtration Device

PlumePort® ActiV® is the first line of defense designed for safe, effective removal of surgical smoke during laparoscopic procedures. ActiV easily attaches to standard trocars and utilizes readily available suction systems providing noiseless operation. With a built-in fluid trap, moisture is effectively managed throughout the entire procedure to ensure uninhibited flow and prevent filter saturation without compromising the integrity of the filtration media.



**INSUFFLATION
SYSTEM**



**PLUMEPORT® SEO
LAPAROSCOPIC**



**PLUMEPORT®
LAPAROSCOPIC**



To learn more about these and other innovative products, call **800-448-6506** or visit **CONMED.com**.



PLUMEPORT® ACTIV®
LAPAROSCOPIC
SMOKE FILTRATION
DEVICE

Features & Benefits

Active Smoke Removal

- Controlled management of smoke throughout the surgical procedure featuring flow control dial with one hand adjustability.

Noiseless Operation

Compact Housing Design

- Contains built-in fluid trap, flow adjustment valve and filtration.

Moisture Management

- Built-in fluid trap to manage moisture throughout entire procedure.

Integrated ULPA Filtration

- Coupled with activated carbon to safely filter peritoneal smoke.

Universal Connectivity

- Easily attaches to standard luer lock fittings and suction systems.

Up to 50% More Flow

- Provides greater flow than passive laparoscopic devices.

Technical Information

Flow Rates

Adjustable flow rates from
3.7 — 12.2 lpm

Based on 5mm trocar with no instruments
with insufflator pressure setting ranging
from 100 to 500

Not made with natural rubber latex.

To learn more about this and other innovative products,
call [1-800-448-6506](tel:1-800-448-6506) or visit CONMED.com.

Ordering Information

	Description	Qty	Catalog Number
NEW	PlumePort® ActiV® with 3/8 in (9.5 mm) x 10 ft (3 m) hose / tubing & 7/8 in (22 mm) filter port adapter	10/Case	PP010CS

Distributed by:
CONMED Corporation
525 French Road, Utica, NY 13502

Customer Service: 1-800-448-6506
International Sales: 1-315-797-8375
Fax: 1-800-438-3051

CONMED.com
info@mail.CONMED.com