





# Specialized Solutions for Proactive Procedures



## A Full Polypectomy Portfolio for Any Budget

Want a simplified polypectomy solution, versatile enough for many clinical scenarios?

**CONMED** has you covered.

CONMED's polypectomy portfolio is designed to be conveniently comprehensive for all your needs. Our commitment to accessibility makes improved outcomes available to every endoscopy suite budget.




Connect with your CONMED representative to trial our polypectomy portfolio today.



# NAKAO Spider-Net®

## Retrieval Device

- Translucent technology\* for optimized visualization
- Unparalleled specimen protection
- Reliable performance without added cost
- Ideal for polyp/tissue retrieval and foreign body extraction



\*When Wet

Product Number	Description	Net Size	Cath Length	Unit /Box
0160MG	NAKAO Spider-Net - EGD	2.5cm x 6cm	160cm	5
00230A	NAKAO Spider-Net		230cm	10
0230AM	NAKAO Spider-Net - Multi-Pack		230cm	20
0350EN	NAKAO Spider-Net Enteroscope		350cm	5

# Optimizer<sup>®</sup> Multi-Chamber Poly Trap

- Four independent, selectable chambers for advanced diagnostic accuracy
- Fast and easy fluid management
- Convenient central vacuum lumen



Product Number	Description	Unit /Box
0160MG	Multi-Chamber Polyp Trap	5
00230A	Multi-Chamber Polyp Trap - Bulk Pack	10


# Singular™ Short Throw Polypectomy Snare

- Short throw handle facilitates one-handed operation and user comfort
- True 2-to-1 pulley system delivers enhanced tactile feel
- Maximum clinical versatility with a full range of loop shapes and sizes

	Loop Stiffness	Product Number	Loop Size	Active Cord Connection	Cath O.D.	Cath Length	Unit /Box
Oval Shape							
 Oval Loop	Firm	000472	Large - 32mm	Olympus	2.3mm	230cm	10
		000474	Med. - 23mm	Olympus			
		000476	Small - 16mm	Olympus			
		000478	X-Small - 11mm	Olympus			
 Crescent Loop	Soft	000472	Large - 32mm	Olympus	2.3mm	230cm	10
		000474	Med. - 23mm	Olympus			
		000476	Small - 16mm	Olympus			
		000478	X-Small - 11mm	Olympus			
Crescent							
 Hexagonal Loop	Firm	000957	Med. - 24mm	Olympus	2.3mm	230cm	10
	Hexagonal						
	Firm	000959	Med. - 25mm	Olympus	2.3mm	230cm	10

# Optimizer<sup>®</sup> Polypectomy Snare

- Gradation to help indicate resection size and progress
- Oval loop is offered in four sizes
- Varying firm and stiffness options




	Loop Stiffness	Product Number	Loop Size	Active Cord Connection	Cath O.D.	Cath Length	Unit /Box
Oval Shape							
Firm		000472	Large - 32mm	Olympus	2.3mm	230cm	10
		000474	Medium - 23mm	Olympus			
		000476	Small - 16mm	Olympus			
		000478	X-Small - 11mm	Olympus			
Soft		000460	Large - 32mm	Olympus	2.3mm	230cm	10
		000462	Medium - 23mm	Olympus			
		000464	Small - 16mm	Olympus			
		000466	X-Small - 11mm	Olympus			

# Orbit-Snare®

## Rotatable Polypectomy Snare


- Rotatable loop enables effective capturing
- Ergonomic handle designed for simple rotation assistance
- Firm wire loop designed for ensnaring flat polyps



Loop Stiffness	Product Number	Loop Size	Active Cord Connection	Cath O.D.	Cath Length	Unit /Box
Oval Shape						
Firm	01-15-230	Small - 15mm	Universal	2.5mm	230cm	10
	01-25-230	Medium - 25mm	Universal			
	01-35-230	Large - 35mm	Universal			

# CONMED Polypectomy Snare

- Most economical line of snares
- Soft and firm loop wire



Loop Stiffness	Product Number	Loop Size	Active Cord Connection	Sheath O.D.	Sheath Length	Unit /Box
----------------	----------------	-----------	------------------------	-------------	---------------	-----------

## Oval Shape

Firm	100600	Small - 16mm	Olympus	2.3mm	230cm	10
	100601	Medium - 23mm	Olympus			
	100602	Large - 32mm	Olympus			
Soft	100603	Small - 16mm	Olympus	2.3mm	230cm	10
	100604	Medium - 23mm	Olympus			
	100605	Large - 32mm	Olympus			





## Precisor<sup>®</sup> EXL<sup>™</sup> Disposable Biopsy Forceps

- Ergonomic, easy-to-use design
- Designed and refined for optimal comfort
- Precise and efficient

Jaw Configuration	Product Number	Active Cord Connection	Jaw O.D.	Shaft Length	Unit /Box
Oval	000852	Olympus	2.3mm	230cm	10
Alligator	000854	Olympus	2.3mm	230cm	10



CONMED Corporation | 11311 Concept Blvd., Largo, FL 33773 | [www.CONMED.com](http://www.CONMED.com)



Service & Support