




Specialized Solutions for Proactive Procedures



A Full Polypectomy Portfolio for Any Budget


Want a simplified polypectomy solution, versatile enough for many clinical scenarios?
CONMED has you covered.



CONMED's polypectomy portfolio is designed to be conveniently comprehensive for all your needs. Our commitment to accessibility makes improved outcomes available to every OR budget.



Connect with your CONMED representative to trial our polypectomy portfolio today.



NAKAO Spider-Net®

Retrieval Device

- Translucent technology* for optimized visualization
- Unparalleled specimen protection
- Reliable performance without added cost
- Ideal for polyp/tissue retrieval and foreign body extraction

*When Wet

Product Number	Description	Net Size	Cath Length	Unit /Box
0160MG	NAKAO Spider-Net - EGD	2.5cm x 6cm	160cm	5
00230A	NAKAO Spider-Net		230cm	10
0230AM	NAKAO Spider-Net - Multi-Pack		230cm	20
0350EN	NAKAO Spider-Net Enteroscope		350cm	5

Optimizer®

Multi-Chamber Polyp Trap

- Four independent, selectable chambers for advanced diagnostic accuracy
- Fast and easy fluid management
- Convenient central vacuum lumen



Product Number	Description	Unit /Box
CAP132	Multi-Chamber Polyp Trap	20
CAP132B	Multi-Chamber Polyp Trap - Bulk Pack	120



Singular™

Short Throw Polypectomy Snare


- Short throw handle facilitates one-handed operation and user comfort
- True 2-to-1 pulley system delivers enhanced tactile feel
- Maximum clinical versatility with a full range of loop shapes and sizes

	Loop Stiffness	Product Number	Loop Size	Active Cord Connection	Cath O.D.	Cath Length	Unit /Box
Oval Shape							
 Oval Loop	Firm	000963	Large - 32mm	Olympus	2.3mm	230cm	10
		000965	Med. - 23mm	Olympus			
		000978	Small - 16mm	Olympus			
		000991	X-Small - 11mm	Olympus			
 Crescent Loop	Soft	000977	Large - 32mm	Olympus	2.3mm	230cm	10
		000983	Med. - 23mm	Olympus			
		000985	Small - 16mm	Olympus			
		000987	X-Small - 11mm	Olympus			
Crescent							
 Hexagonal Loop	Firm	000957	Med. - 24mm	Olympus	2.3mm	230cm	10
		Hexagonal					
	Firm	000959	Med. - 25mm	Olympus	2.3mm	230cm	10

Optimizer®

Polypectomy Snare

- Gradation to help indicate resection size and progress
- Oval loop is offered in four sizes
- Varying firm and stiffness options



Loop Stiffness	Product Number	Loop Size	Active Cord Connection	Cath O.D.	Cath Length	Unit /Box
----------------	----------------	-----------	------------------------	-----------	-------------	-----------


Oval Shape

Firm	000472	Large - 32mm	Olympus	2.3mm	230cm	10
	000474	Med. - 23mm	Olympus			
	000476	Small - 16mm	Olympus			
	000478	X-Small - 11mm	Olympus			
Soft	000460	Large - 32mm	Olympus	2.3mm	230cm	10
	000462	Med. - 23mm	Olympus			
	000464	Small - 16mm	Olympus			
	000466	X-Small - 11mm	Olympus			

Orbit-Snare®

Rotatable Polypectomy Snare


- Rotatable loop enables effective capturing
- Ergonomic handle designed for simple rotation assistance
- Firm wire loop designed for ensnaring flat polyps



Loop Stiffness	Product Number	Loop Size	Active Cord Connection	Cath O.D.	Cath Length	Unit /Box
Oval Shape Firm	01-15-230	Small - 15mm	Universal	2.5mm	230cm	10
	01-25-230	Small - 15mm	Universal			
	01-35-230	Small - 15mm	Universal			

CONMED Polypectomy Snare

- Most economical line of snares
- Soft and firm loop wire



Loop Stiffness	Product Number	Loop Size	Active Cord Connection	Sheath O.D.	Sheath Length	Unit /Box
----------------	----------------	-----------	------------------------	-------------	---------------	-----------

Oval Shape

Firm	100600	16mm	Olympus	2.3mm	230cm	10
	100601	23mm	Olympus			
	100602	32mm	Olympus			
Soft	100603	16mm	Olympus	2.3mm	230cm	10
	100604	23mm	Olympus			
	100605	32mm	Olympus			



Precisor[®] EXL[™]

Disposable Biopsy Forceps

- Ergonomic, easy-to-use design
- Designed and refined for optimal comfort
- Precise and efficient

Jaw Configuration	Product Number	Active Cord Connection	Jaw O.D.	Shaft Length	Unit /Box
Oval	000852	Olympus	2.3mm	230cm	10
Alligator	000854	Olympus	2.3mm	230cm	10



Toll Free: 1-866-4CONMED
International: 727-214-3000

customerexperience@conmed.com
internationalorders@conmed.com

www.CONMED.com

©2022 CONMED Corporation, 891036351 10/22