

# Porta PlumeSafe® 604

## Surgical Smoke Evacuator

Porta PlumeSafe® 604 is an ideal surgical smoke evacuator for veterinarian and non-acute care settings where moderate flow is needed. The unit is portable, easy to use, and features patented technology to track filter life.



**SYSTEM 5000®**  
GENERATOR



**GOLDVAC® SLIM**  
PENCILS



**EZLINK®**  
AUTOMATION ACTIVATION DEVICE



To learn more about these and other innovative products, call **800-448-6506** or visit **CONMED.com**.

## Features & Benefits

### Lightweight & Portable

- Allows for easy transport or tabletop use

### RFID Technology

- Tracks filter life

### LED Display

- Indicates filter life and suction power

### Economical and Easy to Operate

## Technical Information

### Noise Levels (Nominal Values)

@ maximum power < 59 dBA

### Dimensions (H x W x D)

10 in x 6.5 in x 18 in  
25 cm x 17 cm x 46 cm

### Flow Rates

48 cfm (1358 lpm) using 1-1/4 in (32 mm) tubing  
36 cfm (1019 lpm) using 7/8 in (22 mm) tubing

### Weight

With filter installed: 22.0 lbs (9.0 kg)  
Without filter installed: 18.6 lbs (8.0 kg)

### Warranty

One (1) Year Full Warranty



To learn more about this and other innovative products,  
call **1-800-448-6506** or visit **CONMED.com**.

## Ordering Information

	Description	Qty	Catalog Number
<b>NEW</b>	Porta PlumeSafe® 604 (100/120 Volt Systems) filter sold separately	1/Each	PPS604
	EZLink® Sold separately	1/Each	EZLINK01

CONMED Corporation  
525 French Road, Utica, NY 13502

Customer Service: 1-800-448-6506  
International Sales: 1-315-797-8375  
Fax: 1-800-438-3051

CONMED.com  
info@mail.CONMED.com