# CONMED

# **Precisor BRONCHO®**

### Disposable Pulmonary Biopsy Forceps

The Precisor BRONCHO<sup>®</sup> disposable pulmonary biopsy forceps are ideal for obtaining tissue specimens from the airway. A streamlined jaw assembly and high performance handle ensure optimal control and tactile feel, resulting in fast and easy removal of biopsy specimens.

## Features and Benefits

Precise cutting action — large jaw sample capacity — endoscopic channel friendly
Enhanced pushability – reduced wear on biopsy channel
Designed to traverse bends more effectively for easy removal of biopsy specimens
Optimal control and tactile feel – ergonomically designed for comfort

To learn more about these and other innovative products, call **1-866-4CONMED** or visit **www.CONMED.com**.



#### **Design Comparison**

Precisor BRONCHO® vs. Competitive Twin Wire

Jaw Opening Span Comparison






The Precisor BRONCHO® Forceps (left) has a wide jaw span

The Precisor BRONCHO® Forceps (left) has high capacity jaws

#### **Cam Design**



The Precisor BRONCHO® Forceps has internal working components to help maintain biopsy channel integrity.

#### Comparative View from Bronchoscope

Precisor BRONCHO®






Twin Wire Design

Precisor BRONCHO® streamlined jaw assembly provides a clear endoscopic view of the tissue site.

#### **Twin Wire Design**



Accessories with protruding components can cause excessive biopsy channel damage and result in substantial scope repair costs.

## **Ordering Information**

Description	Shaft Length	Jaw O.D	Unit	Catalog Number
Oval Cup (coated)	115cm	1.8mm	10/box	100502
Oval Cup (uncoated)	115cm	1.8mm	10/box	100506
Alligator Cup (coated)	115cm	1.8mm	10/box	100503
Alligator Cup (uncoated)	115cm	1.8mm	10/box	100507

CONMED Corporation 11311 Concept Blvd. Largo, FL 33773 Toll Free: 1-866-4CONMED International: 727-214-3000 www.CONMED.com

# customerexperience@conmed.com internationalorders@conmed.com