



Pressure Infuser Irrigation Pumps

- A full line of single 1-liter and 3-liter systems to meet all of your laparoscopic irrigation pressure source needs
- Quiet and convenient Power Pack Compressor for Infuser Irrigation Pumps




AIRSEAL® SYSTEM
INSUFFLATION MANAGEMENT



UNIVERSAL® PLUS
SUCTION/IRRIGATION SYSTEMS



SPECIMEN BAGS
DISPOSABLE, LAPAROSCOPIC



To learn more about these and other innovative products, call 866-4CONMED or visit www.CONMED.com



1-Liter Pressure Infuser Irrigation Pump

- Takes standard irrigation bag, 500 – 1000ml
- IV pole mountable (3/4" to 1 1/2")
- Locking pressure adjustment knob
- Comes with 12' of gas source tubing to fit connectors below



Dual 1-Liter Pressure Infuser Irrigation Pump

- Takes two bags, 500 – 1000ml, great for irrigators with dual spikes
- IV pole mountable (3/4" to 1 1/2")
- Individual locking pressure adjustment knob for each bag chamber
- Comes with 12' of gas source tubing to fit connectors below



3-Liter Pressure Infuser Irrigation Pump

- Most versatile, takes single bag, 500 – 3000ml
- IV pole mountable (3/4" to 1 1/2")
- Locking pressure adjustment knob
- Comes with 12' of gas source tubing to fit connectors below
- Also can take 2 bags, 500 – 1000ml



Power Pack, Pressure Infuser Irrigation Pump

- AC powered air compressor for 270178, 270123, and 270165 Infuser Irrigation Pumps
- 110V version available
- IV pole mountable (3/4" to 1 1/2")
- UL recognized motor

Ordering Information

Description	Catalog Number
Pressure Infuser Irrigation Pumps	
1-Liter Pressure Infuser Irrigation Pump	270178
Dual 1-Liter Pressure Infuser Irrigation Pump	270123
3-Liter Pressure Infuser Irrigation Pump	270165
Power Pack, Infuser Pump Compressor	270179



Ohmeda



Puritan-Bennett

Gas Source Quick Connect Couplers

Company	Oxygen	Air
Ohmeda	-	CD984A
Puritan-Bennett	CD989	-