





# Expect the Best. We'll do the Rest.

Use state-of-the-art manual instruments  
from the industry leader

CONMED's full line of instruments are made from a solid grooved rod of surgical grade stainless steel for a stronger and more robust feel. Experience the ease and convenience of a single-source supplier for resection instruments.



To learn more about these and other innovative products,  
call **866-4CONMED** or visit **www.CONMED.com**.

# Shutt® Manual Instruments

## FORCEPS

Shutt® has come a long way in setting the standards for arthroscopic instrumentation. As an early pioneer in the field, CONMED produced the first arthroscopy hook punch forceps in 1980.

Shutt® has earned a reputation for precise, dependable instruments. With an unsurpassed commitment to service, Shutt® has met your needs time and time again.

Our instruments feature a barber-tang handle for a sensitive touch and a large thumb ring for long-term comfort. A special coating on the actuator guarantees a smooth cutting action.

Precision-matched, diamond-honed tips assure a clean, close cut.


A heavy-duty shear pin stands up to the toughest tissue. A glass bead finish is easy to clean and gives each instrument a highly-professional look and feel.



### Blunt Tip Forceps

3.4mm dia., 130mm


Straight	1.1001
15° Up	1.1631



### B-Scoop Forceps

3.4mm dia., 130mm

Straight	C6855
30° Left	1.1211
30° Right	1.1221
15° Up	1.1631



### Mini-Tip Forceps

2.75mm dia., 130mm

Straight	2.1003
30° Left	2.1113
30° Right	2.1123
15° Up	2.1633

3.4mm dia., 130mm

15° Up	1.1633
--------	--------



### Shovel Nose Forceps

3.4mm dia., 130mm.


Mini, straight	4.1003
----------------	--------



### Retrograde Forceps

3.4mm dia., 130mm.

Straight, Left Bite	C3007.0
Straight, Right Bite	C3007.1



### 45° Forceps

2.75mm dia., 130mm

Blunt Tip, Right	32.10013
Blunt Tip, Left	32.10014
Scoop Tip, Right	32.10017
Scoop Tip, Left	32.10018

3.4mm dia., 130mm

Blunt Tip, Right	31.10013
Blunt Tip, Left	31.10014
Scoop Tip, Right	31.10017
Scoop Tip, Left	31.10018



### Aggressor® Forceps

3.4mm dia., 130mm

Straight	31.10034
15° Up	31.16334
15° Left	31.16136
15° Right	31.16235



### Mini-Aggressor® Forceps

2.75mm dia., 130mm

Straight	32.10034
15° Up	32.16334

3.4mm dia., 130mm

Straight	31.10033
15° Up	31.16333
30° Left	31.12133
30° Right	31.12233



### Low-Profile Forceps

2.75mm dia., 130mm,  
Radiused Tip

Straight	32.10019
15° Up	32.16319

3.4mm dia., 130mm,  
Radiused Tip

Straight	31.10019
15° Up	31.16319



### Square Trimmer Forceps (SST)

2.75mm dia., 130mm

Blunt Tip, Straight	2.1039
---------------------	--------

Blunt Tip, 15° Up	2.1639
-------------------	--------

3.4mm dia., 130mm, Square

Blunt Tip, Straight	1.1039
---------------------	--------

## SUTURE RETRIEVERS



### Suture Retrieval Forceps

3.4mm dia., 130mm

Straight	16.1018
----------	---------

## SUTURE CUTTERS



### Standard Scissors

2.75mm dia., 130mm

Straight	2.10012
----------	---------

3.4mm dia., 130mm

Straight	1.10012
----------	---------



### Micro-Scissors

2.75mm dia., 130mm

Straight	2.10011
----------	---------

15° Up	2.16311
--------	---------

3.4mm dia., 130mm

Straight	1.10011
----------	---------



### Katana™ Cutters

Flush Cutter	GU1009F
--------------	---------

Suture Cutter	GU1009
---------------	--------

## GRASPERS



### Alligator Grasper

2.75mm dia., 130mm

Straight, w/ratchet	12.1002
---------------------	---------

3.4mm dia., 130mm

Straight, w/ratchet	11.1002
---------------------	---------

4.5mm dia., 130mm

Straight, w/ratchet	18.1002
---------------------	---------



### Grasping Forceps

2.75mm dia., 130mm

Straight, w/ratchet	12.1001
---------------------	---------

3.4mm dia., 130mm

Straight, w/ratchet	11.1001
---------------------	---------

Straight, w/o ratchet	16.1001
-----------------------	---------

4.5mm dia., 130mm

Straight, w/ratchet	18.1001
---------------------	---------

### Custom Grasping Forceps

2.75mm dia., 130mm

Straight, w/o ratchet	17.1001
-----------------------	---------

## SUTURE PASSERS




### Suture Tram®

3.3mm dia., 130mm

Retriever, Straight	C6800
---------------------	-------

Retriever, 15° Upsweft	C6860
------------------------	-------



### Raven™ Suture Passer

3.5mm dia., 138mm

Straight	GU1004
----------	--------

35° Up	GU1008
--------	--------



### Seahawk™ Suture Passer

3.4mm dia., 138mm

Left, 40° Up	GU2000
--------------	--------

Right, 40° Up	GU2001
---------------	--------

RC Left, 30° Down	GU2002
-------------------	--------

RC Right, 30° Down	GU2003
--------------------	--------

## ACCESSORIES

### Instrument Cases

(Shutt Instruments Only)

Basic Set Instrument Case	9899
<i>Accommodates up to fifteen Shutt instruments</i>	

Master Instrument Case	9879
<i>Accommodates up to twenty Shutt instruments</i>	

### Precision Probes

3.5mm dia. tip, 23cm

Mini-Probe, straight, blue handle	21.1001
-----------------------------------	---------

4mm dia. tip, 23cm

Maxi-Probe, straight, green handle	21.1002
------------------------------------	---------


### Forceps Stands

Deluxe Forceps Stand	9898
<i>Accommodates up to fifteen ring-handled instruments</i>	

# Shutt® Specialty Instruments

## SMALL JOINT INSTRUMENTS

Available in 2.0 and 2.7mm diameters and 65 and 90mm lengths, our small joint Shutt® are ideal for arthroscopy in tight spaces. Smaller in size but with the same quality and durability you expect in the Shutt® brand.



### Standard Graspers

2mm dia., 65mm

Straight	15.6001
----------	---------

2.75mm dia., 65mm


Straight	17.6001
----------	---------

2mm dia., 90mm

Straight	15.7001
----------	---------

2.75mm dia., 90mm

Straight	17.7001
----------	---------



### SST Blunt Tip Forceps

2mm dia., 65mm

Straight	5.6009
----------	--------

15° Up	5.60013
--------	---------

15° Down	5.60014
----------	---------

2.75mm dia., 65mm

Straight	2.6009
----------	--------

15° Down	2.60013
----------	---------

2mm dia., 90mm

Straight	5.7009
----------	--------

2.75mm dia., 90mm

Straight	2.7009
----------	--------



### Mosquito Tip Graspers

2mm dia., 65mm

Straight	15.6004
----------	---------

30° Right	15.6005
-----------	---------

2.75mm dia., 65mm

Straight	17.6004
----------	---------

2mm dia., 90mm

Straight	15.7004
----------	---------



### SST Scoop Tip Forceps

2mm dia., 65mm.

Straight	5.6004
----------	--------

15° Up	5.60016
--------	---------


2.75mm dia., 65mm.

15° Down	2.60015
----------	---------

2.75mm dia., 90mm.

Straight	2.7004
----------	--------

15° Down	2.70015
----------	---------



### 45° Forceps

2mm dia., 65mm.

Blunt Tip, Right	35.60013
------------------	----------

Blunt Tip, Left	35.60014
-----------------	----------

2.75mm dia., 65mm.

Blunt Tip, Right	32.60013
------------------	----------

Blunt Tip, Left	32.60014
-----------------	----------

2mm dia., 90mm.

Blunt Tip, Right	35.70013
------------------	----------

Blunt Tip, Left	35.70014
-----------------	----------

## SMALL JOINT ACCESSORIES



### Instrument Case

Small Joint Instrument Case	9896
<i>Accommodates up to ten-ring handled Shutt instruments</i>	

### Probes

65mm Small Joint Probe, 70° tip	26.3003
90mm Small Joint Probe, 70° tip	26.7003
90mm Small Joint Probe, 30° tip	26.7103
Instrument Stabilizer Ball	326.000

### Curettes & Dissectors

2.75mm dia., 65mm

Curette/Dissector	26.3006
-------------------	---------

2.75mm dia., 90mm

Open Curette	22.7005
Dissector	22.7008

## SERIES 210® HIP INSTRUMENTS

Elongated working length (210mm) instruments to enable you to reach into the deep hip joint. Reinforced tubing and welding gives the strength needed to accomplish the toughest hip cases.



### Micro-Scissors

Straight	S210MSC
----------	---------



### Aggressor® Forceps

15° Up	S210AFU
15° Left	S210AFL
15° Right	S210AFR



### Suture Retrieval Forceps

Straight	S210SRF
----------	---------



### Grasping Forceps

Straight w/ratchet	S210GSP
--------------------	---------



### Paradigm® Slim™ Raven™ Suture Passer

15° up	PGU700
--------	--------



### Paradigm® Slim™ Suture Retriever

15° up	PGU701
--------	--------

CONMED Corporation  
11311 Concept Blvd.  
Largo, FL 33773

Toll Free: 1-866-4CONMED  
International: 727-214-3000  
www.CONMED.com

customerexperience@conmed.com  
internationalorders@conmed.com