


SnapEvac™ Surgical Smoke Evacuation Adapter

SnapEvac™ is an ergonomically designed smoke evacuation electrosurgical pencil adapter that transforms most electrosurgical pencils into a surgical smoke evacuation device that captures smoke plume at the point of creation.



SnapEvac™

Shown with ESU pencil.

Shown with monopolar handpiece.

SYSTEM 5000
GENERATOR



GOLDINE®
ESU PENCIL



CLEARVIEW®
SURGICAL SMOKE EVACUATOR



To learn more about these and other innovative products, call **800-448-6506** or visit **CONMED.com**.

Features & Benefits

Snap and Evac

- Just snap in handpiece of your choice

Efficient Workflow

- Eliminates need for OR personnel to hold tubing

Transparent Capture Port

- Improves visibility of the surgical site

Cost-Effective

- Economical option for surgical smoke evacuation

Lightweight Tubing



To learn more about this and other innovative products, call **1-800-448-6506** or visit **CONMED.com**.

Ordering Information

NEW	Description	Qty	Catalog Number
	SnapEvac® Ideal Designed using ABS plastic and polycarbonate capture port with 10 ft (3 m) tubing; 7/8 in (22 mm) cuff. <i>* ESU pencils and handpiece both sold separately. SnapEvac may not fit all ESU pencils and handpieces.</i>	20/Case	SnapEvac20

CONMED Corporation
525 French Road, Utica, NY 13502

Customer Service: 1-800-448-6506
International Sales: 1-315-797-8375
Fax: 1-800-438-3051

CONMED.com
info@mail.CONMED.com