



Specimen Bags

A COMPREHENSIVE OFFERING



ANCHOR* TISSUE RETRIEVAL SYSTEM™

“Ask for the purple one”™

12 RIP-STOP NYLON OPTIONS

Simple

- Designed for ease of use
- Opens and unrolls automatically when deployed

Strong

- Rip-stop nylon for exceptional strength and quality when compared to polymer
- Provides strong and secure seal for specimen retrieval

Innovative

- Unique deployment mechanism offers safety and security
- Proprietary bag material resists breaks and leaks
- Bag shapes engineered to accommodate multiple surgical specialties
- Redeployable



CONMED SPECIMEN BAGS

Simple and effective specimen retrieval

DETACHABLE LAPAROSCOPIC SPECIMEN BAGS

- Disposable specimen bags are supplied with an introducer and bag for ease of insertion through a 5, 10, or 12mm trocar cannula
- Insert indicator ensures bag will be in proper position when deployed
- Simple, ergonomic, one-handed operation

INDEPENDENT MEMORY WIRE BAGS

- Disposable specimen bags with memory wire in a dispensing tube for insertion through a standard 10 or 12mm cannula





Connect with your local representative today or call 1-866-4CONMED.

ORDERING & PRODUCT INFORMATION

Trocar Size	Mouth Diameter	Bag Length	Volume	Qty	Catalog Number
Anchor* - Rip-Stop Nylon Retrieval System					
8mm	2in (5.1cm)	5.5in (14cm)	125mL	5/Box	TRS-ROBO-8
8mm	2.2in (5.6cm)	5.5in (14cm)	175mL	5/Box	TRS-ROBO-8L
10mm	2.2in (5.6cm)	5.5in (14cm)	175mL	5/Box	TRS100SB2
10mm	3in (7.6cm)	5.5in (14cm)	250mL	5/Box	TRS100SB2L
12mm	3in (7.6cm)	5.5in (14cm)	250mL	5/Box	TRS-ROBO-12
15mm	5.5in (14cm)	9.9in (25.4cm)	1200mL	3/Box	TRS-VATS-15
15mm	5.5in (14cm)	9.9in (25.4cm)	1200mL	3/Box	TRS175SB2
15mm	6in (15.2cm)	10.4in (26.5cm)	1800mL	3/Box	TRS190SB2
25mm	6.5in (16.5cm)	15in (38.1cm)	6000mL	3/Box	TRS-TV-25
Anchor* - Rip-Stop Nylon Bag Alone					
10mm	2.2in (5.6cm)	5.5in (14cm)	175mL	10/Box	TRS070
10mm	3.5in (8.9cm)	5.5in (14cm)	250mL	10/Box	TRS080
15mm	5.5in (14cm)	9.9in (25.4cm)	1200mL	5/Box	TRS090
Nylon - Attached Memory Wire Bag					
10mm	3in (7.6cm)	6in (15.2cm)	240mL	5/Box	SB1036
12mm	7.5in (19.1cm)	9in (22.9cm)	1800mL	5/Box	SB1079
Polymer - Attached Memory Wire Bag					
5mm	3in (7.6cm)	4in (10.2cm)	120mL	10/Box	SB534
10mm	3in (7.6cm)	6in (15.2cm)	330mL	5/Box	SB936
10mm	5in (12.7cm)	7in (17.7cm)	750mL	5/Box	SB957
12mm	7.5in (19.1cm)	9in (22.9cm)	1800mL	5/Box	SB979
Polymer - Independent Memory Wire Bag					
10mm	3in (7.6cm)	6in (15.2cm)	200mL	10/Box	SB836
10mm	5in (12.7cm)	7in (17.7cm)	750mL	10/Box	SB857

* Anchor is a registered trademark of Addison Anchor Product Co.

CONMED Corporation
11311 Concept Blvd.
Largo, FL 33773

Toll Free: 1-866-4CONMED
International: 727-214-3000
www.CONMED.com

customerexperience@conmed.com
internationalorders@conmed.com