



## Stat 2® & Stat 2® Pumpette I.V. Gravity Flow Controllers

Offering convenience and safety to both patients and healthcare providers, Stat 2® provides control of I.V. fluid and non-critical medications when electronic infusion is not medically necessary.



**PADPRO®**  
ADULT ELECTRODES



**PULSE OXIMETRY**  
COMPATIBLE SENSORS



**VENI-GARD® I.V.**  
STABILIZATION DRESSING



To learn more about these and other innovative products, call **800-448-6506** or visit **ConMed.com**.



STAT 2®  
20 DROP



STAT 2®  
60 DROP



STAT 2®  
12" EXTENSION SET



# Stat 2® I.V. Gravity Flow Control Device

Available in full primary sets for use in various settings such as Acute Care, Sub-Acute Care and Home Care, the Stat 2® provides excellent control of I.V. fluid and non-critical medications.

## Key features include:

- One handed operation
- Security window with snap lock minimizing accidental manipulation of dial during therapy.
- Ideal for routine hydration and antibiotic therapies
- Needle-Free versions

## Stat 2® & Stat 2® Pumpette Control Devices are available in:

- 20 and 60-drop primary sets with a 15-micron filter within the universal drip chamber, back-check valve, a distal luer/slip lock connector and 84" length.
- Either two split septum "Y" sites, or two swabbable valve "Y" sites.
- Extension sets of 12" with one "Y" site, either split septum or swabbable valve.

# Stat 2® Pumpette Compensating I.V. Gravity Flow Control Device

Easy and intuitive to use, the Stat 2® Pumpette's unique design compensates for changes in both I.V. bag volume, as well as in bag height.

**Key features include:**

- Pumpette controller can be set at eye level for convenient clinical verification and adjustment
- Security window minimizes accidental manipulation of dial face during therapy
- Ideal I.V. delivery system for routine hydration and antibiotic therapies
- Needle-Free versions



STAT 2® PUMPETTE  
20 DROP



STAT 2® PUMPETTE  
60 DROP



STAT 2® PUMPETTE  
12" EXTENSION SET

I.V. GRAVITY  
FLOW  
CONTROLLERS

*All I.V. controllers contain slide clamp, distal luer / slip lock connection.*  
*All primary sets include a back check valve and universal spike with a 15 micron filter.*  
*These devices are not made with natural rubber latex.*



525 French Road, Utica, NY 13502

I.V. GRAVITY  
FLOW  
CONTROLLERS

Ordering Information

Description	Quantity		Catalog Number
<b>Stat 2® I.V. Gravity Flow Controllers</b>	<b>Pouch/Pkg</b>	<b>Case</b>	
12" length extension set, controller, 1 split-septum "Y" site, distal luer / slip lock connector	1	60	S2-12
20 drop/mL, 84" primary set controller, 2 split-septum "Y" sites, universal spike, 15-micron filter, back-check valve, distal luer / slip lock connector	1	48	S2V-20
60 drop/mL, 84" primary set controller, 2 split-septum "Y" sites, universal spike, 15-micron filter, back-check valve, distal luer / slip lock connector	1	48	S2V-60
12" length extension set, controller, 1 swabbable valve "Y" site, distal luer / slip lock connector	1	60	S2-12 N
20 drop/mL, 84" primary set controller, 2 swabbable valve "Y" sites, universal spike, 15-micron filter, back-check valve, distal luer / slip lock connector	1	48	S2V-20 N
60 drop/mL, 84" primary set controller, 2 swabbable valve "Y" sites, universal spike, 15-micron filter, back-check valve, distal luer / slip lock connector	1	48	S2V-60 N
<b>Stat 2® Pumpette I.V. Gravity Flow Controllers</b>	<b>Pouch/Pkg</b>	<b>Case</b>	
12" length extension set, pumpette controller, 1 split-septum "Y" site, distal luer / slip lock connector	1	60	P-S2-12
20 drop/mL, 84" primary set, pumpette controller, 2 split-septum "Y" sites, universal spike, 15-micron filter, back-check valve, distal luer / slip lock connector	1	48	P-S2V-20
60 drop/mL, 84" primary set, pumpette controller, 2 split-septum "Y" sites, universal spike, 15-micron filter, back-check valve, distal luer / slip lock connector	1	48	P-S2V-60
12" length extension set, pumpette controller, 1 swabbable valve "Y" site, distal luer / slip lock connector	1	60	P-S2-12 N
20 drop/mL, 84" primary set, pumpette controller, 2 swabbable valve "Y" sites, universal spike, 15-micron filter, back-check valve, distal luer / slip lock connector	1	48	P-S2V-20 N
60 drop/mL, 84" primary set, pumpette controller, 2 swabbable valve "Y" sites, universal spike, 15-micron filter, back-check valve, distal luer / slip lock connector	1	48	P-S2V-60 N
<b>Accessories</b>			
Infant size 5.25" x 2.25" arm support with covers	60 supports, 90 covers		810-L
Pediatric size 7" x 3" arm support with covers	60 supports, 90 covers		910-L
Adult size 10.5" x 4" arm support with covers	80 supports, 80 covers		1010-L
Adult size arm supports only	80 supports		1011-L
Adult size arm supports covers	250 covers		1012

Customer Service: 1-800-448-6506  
International Sales: 1-315-797-8375  
Fax: 1-800-438-3051

ConMed.com  
info@mail.conmed.com