


MicroNeedle

MicroNeedle

ConMed's MicroNeedles are designed to fit most standard electrosurgical pencils that accommodate an electrode with a $\frac{3}{32}$ inch (2.4 mm) or $\frac{5}{32}$ inch (4 mm) shaft.

- Molded with rings to help grip electrodes
- Ultra fine tip
- Heat-resistant tungsten
- Use at low power settings



Ordering Information

Catalog Number	Product Description	Packaging
138504-10	5 cm, 45° MicroNeedle	10/case 1/pkg
138504	5 cm, 45° MicroNeedle	40/case 1/pkg
138505-10	3 cm MicroNeedle	10/case 1/pkg
138505	3 cm MicroNeedle	40/case 1/pkg
138506-10	4 cm MicroNeedle	10/case 1/pkg
138506	4 cm MicroNeedle	40/case 1/pkg

Please contact your local ConMed Territory Manager for information on our complete product line and accredited training programs.

Manufactured by: **Unimed**
Distributed by: ConMed
Made in USA

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road
Utica, NY 13502
(315) 797-8375
Customer Service: 1-800-448-6506
FAX: 1-800-438-3051
International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235
email: info@mail.conmed.com
www.conmed.com

