

ULPA

Surgical Smoke Evacuator Replacement Filter

For the Stryker NEPTUNE 2 and NEPTUNE 3



PLUMEPEN® ULTRA
PENCILS



PLUMEPEN® ELITE
PENCILS



GOLDVAC® SLIM
PENCILS



To learn more about these and other innovative products, call **800-448-6506** or visit **CONMED.com**.



ULPA
SURGICAL SMOKE
EVACUATOR
REPLACEMENT
FILTER

Features & Benefits

Compatibility

- Meets same performance criteria of original equipment manufacturer

Filter Life

- 80 hours

ULPA Filtration

- Safely filters out particles .1 micron and larger produced during surgical procedures

Odor Management

- Integrated activated carbon may help to reduce noxious odors

Multi-Port Design

- Multi-port reduces or eliminates the need for adapters when using various accessories

Technical Information

Filter Performance

ULPA Filter .1 to .2 particles captures at 99.999% efficiency

Dimensions (H x W x D)

7.0 in x 7.4 in x 7.9 in
17.78 cm x 18.796 cm x 22.86 cm

Filter Ports

1/4 in x 3/8 in x 7.9 in
6.4 mm x 9.5 mm x 22 mm

All other performance criteria not listed here are equivalent to original equipment manufacturer part number 0702-040-000.

To learn more about this and other innovative products, call **1-800-448-6506** or visit **CONMED.com**.

Ordering Information

	Description	Qty	Catalog Number
NEW	ULPA Surgical Smoke Evacuator Replacement Filter for the Stryker NEPTUNE 2 and NEPTUNE 3	1/Each	VSNEP01
	ULPA Surgical Smoke Evacuator Replacement Filter for the Stryker NEPTUNE 2 and NEPTUNE 3	4/Each	VSNEP04

Distributed by:
CONMED Corporation
525 French Road, Utica, NY 13502

Customer Service: 1-800-448-6506
International Sales: 1-315-797-8375
Fax: 1-800-438-3051

CONMED.com
info@mail.CONMED.com