





VCare[®] Plus

Cervical Cup Sizing



Size	Catalog Number	Cervical Cup Diameter	Cervical Cup Volume	Cervical Cup Nominal Depth
Small (S) 	60-6085-200A	32 mm (1.26 in)	9.7 cm ³ (0.59 in ³)	0.660 in
Medium (M) 	60-6085-201A	34 mm (1.34 in)	14.3 cm ³ (0.87 in ³)	0.876 in
Large (L) 	60-6085-202A	37 mm (1.46 in)	20.7 cm ³ (1.26 in ³)	1.053 in
Extra Large (XL) 	60-6085-203A	40 mm (1.57 in)	26.7 cm ³ (1.63 in ³)	1.230 in