



## VCare<sup>®</sup> Plus and VCare<sup>®</sup> DX Plus

VCare<sup>®</sup> Plus helps surgeons deliver clear control of the anatomy using effective uterine manipulation while assisting with the identification of key surgical landmarks for a safe, efficient, and repeatable procedure.



**AIRSEAL<sup>®</sup>**  
ACCESS MANAGEMENT SYSTEM



**ANCHOR<sup>\*</sup>**  
TISSUE RETRIEVAL SYSTEM™



**VISICLEAR<sup>®</sup>**  
SURGICAL SMOKE EVACUATOR



To learn more about these and other innovative products, call **866-4CONMED** or visit **www.CONMED.com**.

\*Anchor is a trademark of the Anchor Products Company, Addison, Illinois.


## ➤ VCare® uterine manipulators help you deliver a better surgical experience.

- Provides superior exposure and access to pelvic cavity
- Delineates vaginal fornices for easy identification of colpotomy plane
- Maintains pneumoperitoneum throughout the procedure



**VCare® DX Plus provides excellent control of the uterus without distention of vaginal fornices. Ideal for:**

- Tubal Ligation
- Supracervical Hysterectomy
- Diagnostic Laparoscopy
- Tubal Occlusion



To learn more about this and other innovative products, call **866-4CONMED** or visit [www.CONMED.com](http://www.CONMED.com).

## Ordering Information

Description	Qty	Catalog Number
<b>VCare® Plus</b>		
Small (32mm) Cervical Cup	8/Case	60-6085-200A
Medium (34mm) Cervical Cup	8/Case	60-6085-201A
Large (37mm) Cervical Cup	8/Case	60-6085-202A
X-Large (40mm) Cervical Cup	8/Case	60-6085-203A
<b>VCare® DX Plus</b>		
VCare® DX Plus Manipulator	8/Case	60-6080-000A