

# VisiClear<sup>®</sup>

## Surgical Smoke Evacuator

VisiClear<sup>®</sup> is our most advanced acute care offering for surgical smoke evacuation. Combining Safeport Technology<sup>™</sup>, an occlusion warning, filter life tracking, and one-touch programmable procedure modes for open tubing, electrosurgical pencil, and laparoscopic procedures, VisiClear effectively captures and filters surgical smoke, removes odors, particulates and other potentially hazardous by-products.



**SYSTEM 5000**  
GENERATOR



**PLUMEPEN<sup>®</sup> ULTRA**  
PENCILS



**EZLINK<sup>®</sup>**  
AUTOMATION ACTIVATION DEVICE



**VIROSAFE<sup>®</sup> FLUID TRAP**  
FILTER



To learn more about these and other innovative products, call **866-4CONMED** or visit **www.CONMED.com**.



**BUFFALO FILTER®**

**VISICLEAR®  
SURGICAL SMOKE  
EVACUATOR**

## Features & Benefits

### Procedure Operating Modes

- One of three operating modes (Laparoscopic, Open Tubing, and Electrosurgical Pencil Mode) can easily be selected, automatically customizing filter life based on surgical procedure.

### Advanced Sound Control Design

- Whisper Technology™ ensures ultra-quiet operation with unprecedented sound quality.

### Sleek Shelf Design

- Ergonomically designed to accommodate stacking.

### Occlusion Warning

- Safety feature that alerts medical personnel when an occlusion has occurred.

### Intuitive Color Touch Screen

- Easy to use, requires minimal training, and can be operated with gloved hand.

### Unique Sliding Port Closure

- SafePort Technology™ ensures reduced risk of staff coming in contact with a contaminated port.

### Technical Information

#### Noise Levels (Nominal Values)

@ minimum power < 43 dBA

@ maximum power < 55 dBA

#### Dimensions (H x W x D)

7.1 in x 13.9 in x 16.1 in

18.0 cm x 35.4 cm x 40.9 cm

#### Flow Rates

30 cfm (839 lpm) using 7/8 in (22 mm) tubing

#### Weight

With filter installed: 17.5 lbs (7.9 kg)

Without filter installed: 14.5 lbs (6.6 kg)

#### Static Suction

150 mmHg vacuum

#### Certifications

UL, cUL, CE

#### Warranty

Three (3) Year Full Warranty

To learn more about this and other innovative products, call **866-4CONMED** or visit [www.CONMED.com](http://www.CONMED.com).

## Ordering Information

Description	Qty	Catalog Number
VisiClear® (100/120 Volt System) Includes one filter	1/Each	VC120
VisiClear® (200/240 Volt System) Includes one filter	1/Each	VC220
VisiClear® (100/120 Volt System) Does not include filter	1/Each	VC120 NF

CONMED Corporation  
11311 Concept Blvd, Largo, FL 33773

866-4CONMED  
International: 727-214-3000

[www.CONMED.com](http://www.CONMED.com)  
[Customerexperience@CONMED.com](mailto:Customerexperience@CONMED.com)