



## Click-Tip™ Injection Needles

Designed for controlled needle deployment and injection, CONMED offers a full line of sclerotherapy needles to meet your needs.

CONMED's Click-Tip™ needle family includes 25, 22 and 19 gauge needles for accommodating a variety of clinical situations. The ergonomic, nurse-designed handle is easy to use for maximum comfort.



## Features and Benefits

### Top Down Injection

Allows for easy injection of solutions while the parallel port enhances handling and injection ergonomics.

### Multiple Gauges and Lengths

CONMED has a full line of needles to meet your needs and accommodate a variety of procedure types.

### Color Indication on "Unlock Button"

Clearly shows how to unlock needle in a low-light setting.

### Needle Control System Designed for Safety

Safety Tab decreases the risk of scope damage and limits risk of sharps injury while the locking feature helps minimize needle bounce.

### 15° Beveled Stainless Steel Needle

Provides maximum piercing capabilities for easy insertion.



**To arrange a demonstration, call 800-448-6506.  
To locate your local CONMED Sales Representative  
or visit [CONMED.com](http://CONMED.com) for more information.**

To learn more about this and other innovative products, call  
**1-800-448-6506** or visit **CONMED.com** for more information.

## Ordering Information

Description	Needle Length	Catheter O.D.	Catheter Length	Qty	Catalog Number
<b>CLICK-TIP™ Injection Needle</b>					
Click-Tip™ Injection Needle, 19 Gauge	4mm	2.3mm	180cm	10/box	02-23-180
Click-Tip™ Injection Needle, 19 Gauge	4mm	2.3mm	230cm	10/box	02-23-230
Click-Tip™ Injection Needle, 19 Gauge	6mm	2.3mm	180cm	10/box	04-23-180
Click-Tip™ Injection Needle, 19 Gauge	6mm	2.3mm	230cm	10/box	04-23-230
Click-Tip™ Injection Needle, 22 Gauge	4mm	1.9mm	180cm	10/box	01-19-180
Click-Tip™ Injection Needle, 22 Gauge	4mm	1.9mm	230cm	10/box	01-19-230
Click-Tip™ Injection Needle, 22 Gauge	6mm	1.9mm	180cm	10/box	03-19-180
Click-Tip™ Injection Needle, 22 Gauge	6mm	1.9mm	230cm	10/box	03-19-230
Click-Tip™ Injection Needle, 25 Gauge	4mm	1.9mm	180cm	10/box	05-19-180
Click-Tip™ Injection Needle, 25 Gauge	4mm	1.9mm	230cm	10/box	05-19-230
Click-Tip™ Injection Needle, 25 Gauge	6mm	1.9mm	180cm	10/box	06-19-180
Click-Tip™ Injection Needle, 25 Gauge	6mm	1.9mm	230cm	10/box	06-19-230

Customer Service: 1-800-448-6506  
International Sales: 1-315-797-8375  
Fax: 1-800-438-3051

CONMED.com  
info@mail.CONMED.com