




NEW

## UltrAblator™ Bipolar Probe Series

The NEW UltrAblator™ Bipolar Probe Series utilizes your existing monopolar ESU to provide all the clinical benefits of bipolar ablation without any new capital investment. These unique probes allow you to achieve stable ablation performance with low power output designed to maintain low temperatures, helping keep your patients and equipment safe.



**24K® PUMP**  
FLUID MANAGEMENT SYSTEM



**EDGE™**  
BIPOLAR ARTHROSCOPIC RF SYSTEM



**TRUE HD 3MOS**  
CAMERA SYSTEM



To learn more about these and other innovative products, call **800-237-0169** or visit **CONMED.com**.



90° probe  
for Shoulder Procedures

ULTRABLATOR™  
BIPOLAR PROBE  
SERIES

## Introducing bipolar ablation without dedicated capital

With the UltrAblator™ Bipolar Probe Series, you can finally have clinically-effective bipolar ablation without having to invest in new capital equipment or endure high per procedure costs:

### Clinical Benefits

#### Effective Bipolar Ablation

- Quick, controlled plasma-based ablation and coagulation of tissue in arthroscopic environments
- Optional aspiration designed to improve visibility\*
- Uses your ESU's existing foot pedal\*\*

#### Low Power Output, Low Temperature

- UltrAblator™ 3.3mm, 90° probe generates lower heat than existing systems in the market<sup>1</sup>
- Helps mitigate delays caused by flushing the joint or waiting for fluid to cool down
- Ideal for surgeons concerned with high-joint temperatures

#### Portfolio to Meet Wide Range of Procedures

- Shoulder: 3.3mm, 90° – non-aspirating & aspirating probes
- Knee: 3.3mm, 50° – non-aspirating & aspirating probes
- Small Joint: 2.4mm, 50° – non-aspirating probes, 1.8mm, 60° – non-aspirating probes

<sup>1</sup> Data on file. \*Suction only available on 3.3mm electrodes. \*\*UltrAblator Bipolar Probe Series require a universal adapter and an ESU with monopolar foot control capabilities.

# What if lowering per procedure costs was as easy as switching probes?

## Benefits for OR Staff, Materials Managers and Purchasing

### No Dedicated Capital Equipment Required

- Works with most standard monopolar generators\*\* eliminating new capital investments
- Less equipment to purchase, service, store and maintain
- Saves space in your OR
- Eliminates need to train staff on dedicated bipolar generator
- Makes incorporating bipolar technology simple and affordable

### Lower Per Procedure Costs

- Lower probe costs than many expensive high-performance probes on the market
- No patient grounding pad required

*UltraAblator™ Bipolar Probe Series use a monopolar ESU to provide the clinical benefits of bipolar ablation.*

**Looking for an Aggressive Ablation Solution?  
Check Out CONMED's Edge™ Bipolar Arthroscopic  
RF System at [CONMED.com](http://CONMED.com).**

\*\*UltraAblator Bipolar Probe Series require a universal adapter and an ESU with monopolar foot control capabilities.

ULTRABLATOR™  
BIPOLAR PROBE  
SERIES



COMPATIBLE WITH  
SYSTEM 2450



50° probe  
for Knee and Small  
Joint Procedures



ULTRABLATOR™  
BIPOLAR PROBE  
SERIES

*To take the next step toward lowering the costs of bipolar technology, call **800.237.0169** to locate your local CONMED Sales Representative or visit **CONMED.com** for more information.*

## Ordering Information

Description	Diameter	Qty.	Catalog No.
<b>UltrAblator™ Bipolar Probe Series</b>			
UltrAblator™ Bipolar Ablator, 90° Aspirating Probe	3.3mm	1	C3390A
UltrAblator™ Bipolar Ablator, 90° Non-Aspirating Probe	3.3mm	1	C3390NA
UltrAblator™ Bipolar Ablator, 50° Aspirating Probe	3.3mm	1	C3350A
UltrAblator™ Bipolar Ablator, 50° Non-Aspirating Probe	3.3mm	1	C3350NA
UltrAblator™ Bipolar Ablator, 55° Non-Aspirating Probe	2.4mm	1	C2455NA
UltrAblator™ Bipolar Ablator, 60° Non-Aspirating Probe	1.8mm	1	C1860NA
<b>Electrosurgical Unit (ESU)</b>			
System 2450™ ESU, 100V	–	1	60-2450-100
System 2450™ ESU, 120V	–	1	60-2450-120
System 2450™ ESU, 220V	–	1	60-2450-220
System 2450™ ESU, 230-240V	–	1	60-2450-230
<b>Accessories</b>			
Universal® Adapter	–	1	60-6790-001

CONMED Corporation  
525 French Road, Utica, NY 13502

Local 727-392-6464  
Toll Free 1-800-237-0169

CONMED.com  
customer\_service@conmed.com