




ADVANCED ENDOSCOPIC TECHNOLOGIES

Gastroenterology and Pulmonology Product Catalog




DURACLIP®
HEMOSTASIS CLIP



OPTIMIZER®
POLYPECTOMY SNARE



ELIMINATOR® PET
BALLOON DILATOR



BEAMER AVEO™
ELECTROSURGICAL PLATFORM



To order these and other
innovative products,
call **866-4CONMED** or
visit www.CONMED.com.



TABLE OF CONTENTS

ORDERING INFORMATION

Ordering Information	1
----------------------	---

BILIARY

Papillotomes

TruPass® SW Biliary Papillotome	2
Apollo 3AC® Advanced Cannulation Triple Lumen Papillotomes	3
Apollo 3® Tapered Tip Triple Lumen Papillotomes	3
CompleteControl™ Sphincterotome	4

Cannulae

ProForma® HF 4.5 Cannula	4
--------------------------	---

Guidewires

FXWire® Advanced Measurement Guidewire	5
XWire® Guidewire	6

Dilation Balloons

Eliminator® PET Biliary Balloon Dilators	7
--	---

Stone Removal Balloons

FIT® Biliary Stone Balloon	8
DuraGlide™ Triple Lumen Stone Removal Balloons	8

ENDOSCOPIC ULTRASOUND

Fine Needle Aspiration

ClearView™ Endoscopic Ultrasound Needle	9
---	---

DILATION

PET Balloon Dilators

Eliminator® PET Balloon Dilators	10
CONMED Balloon Inflation Device	11
Multi-Stage Inflation Pressures	11

Polyvinyl Dilators

CleanGuide™ Disposable OTW Esophageal Dilators	12
CleanGuide™ Disposable Spring Tip Guidewire	13

HEMOSTASIS

Hemostasis Clips

DuraClip® Hemostasis Clip	14
---------------------------	----

Hemostasis Probes

BiCap® Superconductor™ Multi-electrode Bipolar Probes	15
Beamer AVEO™ Argon Probes	16

Disposable Injection Needles

Click-Tip™ Injection Needle	17
Flexitip™ Disposable Sclerotherapy Needles	17

Ligation Devices

Stiegmann-Goff™ Endoscopic Ligator Kit	18
Stiegmann-Goff™ Bandito™ Endoscopic Ligator	19

BIOPSY

Disposable Biopsy Forceps

Precisor® EXL™ Coated Disposable Biopsy Forceps	20
Optibite® Disposable Biopsy Forceps	20
Pediatric GI	21
Precisor® Hot Disposable Biopsy Forceps	21
Precisor® Broncho® Disposable Biopsy Forceps	21

Disposable GI Cytology Brushes

CONMED Cytology Brushes	22
Cytology Brushes	22

POLYPECTOMY

Disposable Snares

Beamer® Argon Snare Probe	23
Singular® Polypectomy Snares	24
Optimizer® Polypectomy Snares	25
Orbit-Snare® Rotatable Polypectomy Snares	26
CONMED Polypectomy Snares	26
Active Cord	26

Specimen Devices

Optimizer® Multi-Chamber Polyp Trap	28
-------------------------------------	----

INFECTION PREVENTION

Biopsy Valves / Introducer

CleanGuard® Biopsy Valves	29
Biopsy Channel Accessory Introducer	29

Bite Blocks

Scope Saver® Bite Blocks	29
CONMED Bite Blocks	29

Endoscope Valve-Tracking and Pre-Cleaning System

CleanCage™ Endoscope Valve-Tracking System	30
CleanTray™ Endoscope Pre-Cleaner	30

CleanGuard™ Disposable Endoscope Valves

CleanGuard™ Disposable Endoscope Valves	31
---	----

Cleaning Brushes

Channel Master™ Endoscope Cleaning Brushes	32
Whistle® Endoscope Cleaning Brush	32

PULMONARY

Transbronchial Aspiration Needles

WANG® Cytology Needles	34
WANG® Histology Needles	34

Biopsy Forceps

Precisor Broncho® Disposable Biopsy Forceps	35
---	----

Brushes

Cytology Brushes	36
Needle Tip and Peripheral Cytology Brushes	37
Microbiology Brush	37

Cleaning Brushes

Channel Master™ Endoscope Cleaning Brush	38
Reusable Cleaning Brush	38

GENERATORS & ARGON PROBES

Generators & Argon Probes

Beamer AVEO™	39
Beamer AVEO™ Argon & SideFire Probes	39
Beamer® System CE600	39
Beamer® Argon Probe	40
Beamer® Argon SideFire Probe	40
Beamer® Argon Snare Probe	40
BiCap® III Electrosurgical Generator	41
AVEO Pump	41
Coaxial Bipolar Adapters & Active Cord	42
ABCFlex® Argon Probes	42

LATEX PRODUCTS 43

CATALOG NUMBER INDEX 44

To Place Your Order

Please provide a purchase order number, billing and shipping instructions, product catalog number, description and pricing to:

Customer Service if calling from within the U.S.

Telephone: (866) 4CONMED
Fax: (800) 438-3051
E-mail: customerservice@CONMED.com

International Customer Service

Telephone: +1(727) 214-3000
Fax: +1(315) 735-6235

Terms and Conditions

All prices and specifications are subject to change without notice. All prices are F.O.B. Utica, New York and Lithia Springs, Georgia, (F.O.B. Drop Shipment Points where applicable). Purchase orders should be addressed to CONMED Corporation, 525 French Road, Utica, NY 13502.

CREDIT TERMS: Net 30 days.

Return Goods Policy

All CONMED products, except products that are obsolete, discontinued, or damaged and specialty products, may be returned within 60 days after date of purchase for credit subject to the following:

1. A copy of the invoice number and date of purchase and a Return Authorization Number must accompany the returned goods with the reason for the return stated in writing. Contact CONMED Customer Service for the Return Authorization Number at 1-866-4CONMED.
2. All returned goods would be subject to a minimum of 15% restocking charge. The exact amount of credit to be given, if any, will be determined by a Quality Assurance inspection. Sterile packaged products may be returned for credit only if returned in the original, unopened shelf package.
3. If the return is due to our error in shipping, full credit will be given.
4. All returns should be addressed as indicated by Return Authorization.

Scheduled Shipments

To help maximize the efficiency of your ordering process and inventory control, we offer scheduled shipments.

The advantages include:

- regular delivery schedule
- one purchase order per year
- favorable pricing
- price protection
- free freight for Federal Express Economy Delivery

Quality Policy

At CONMED our passion and responsibility is for the success of our customers and the patients for whom they care.

We empower our customers to succeed by our commitment to integrity, quality, responsiveness and a relentless focus to deliver accessible CONMED solutions. Our talented employees and culture are the foundation of our success.

Papillotomes



TruPass® Short Wire Triple Lumen Papillotome

Sterile, Disposable

Description	Catheter Length	Tip O.D.	Catheter O.D.	Cutting Wire Length	Cutting Wire Type	Guidewire Compatibility	QTY	Catalog Number
TruPass® Triple Lumen Papillotome								
TruPass®	190cm	4.5 F	7-4.5 F	20mm	Monofilament	035" / .89mm	1/box	T2043F
TruPass®	190cm	4.5 F	7-4.5 F	30mm	Monofilament	035" / .89mm	1/box	T3043F
TruPass® SW Triple Lumen Papillotome with Guidewires								
TruPass®with FXWire®	190cm	4.5 F	7-4.5 F	20mm	Monofilament	035" / .89mm	1/box	P2043F
TruPass®with FXWire®	190cm	4.5 F	7-4.5 F	30mm	Monofilament	035" / .89mm	1/box	P3043F

- Guaranteed Cutting Wire Orientation between 11:00 - 1:00
- 4.5 F .035" Compatible Triple Lumen Tapered Tip
- Smooth, Atraumatic Radiopaque Tip
- Easily identifiable green guidewire exit port
- Exit Marker Band
- Cutting Wire Markers

Papillotomes




Apollo 3AC® Advanced Cannulation Triple Lumen Papillotomes

Sterile, Disposable

Catalog Number	Type	Tip		Cutting Wire		Catheter		Guidewire	
		Length	O.D.	Length	Type	Length	O.D.	Compatibility	Unit
7104AC	Short Nose	5mm	4.5 F	20mm	Monofilament	190cm	7-4.5 F	.035"/.89mm	1/box
7105AC	Short Nose	5mm	4.5 F	30mm	Monofilament	190cm	7-4.5 F	.035"/.89mm	1/box

- Guaranteed Cutting Wire Orientation
- 4.5 F .035" Compatible Triple Lumen Tapered Tip
- Color Coded Guidewire Port
- Exit Marker Band
- Cutting Wire Markers




Apollo 3® Tapered Tip Lumen Papillotomes

Sterile, Disposable

Catalog Number	Type	O.D.	Cutting Wire		Catheter		Guidewire	
			Length	Type	Length	O.D.	Compatibility	Unit
007104	Short Nose	5.0 F	20mm	Monofilament	190cm	7-5.0 F	.035"/.89mm	1/box
007105	Short Nose	5.0 F	30mm	Monofilament	190cm	7-5.0 F	.035"/.89mm	1/box

Papillotomes



CompleteControl™ Sphincterotome

Catalog Number	Description	Catheter Length	Insulated Cutting Wire Length	Cutting Wire Type	QTY
SLW20-5096-5	CompleteControl™ Sphincterotome, Long Wire	180cm	20mm	Monofilament	5/box
SLW30-5096-5	CompleteControl™ Sphincterotome, Long Wire	180cm	30mm	Monofilament	5/box
SSW20-50912-5	CompleteControl™ Sphincterotome, Short Wire	180cm	20mm	Monofilament	5/box
SSW30-50912-5	CompleteControl™ Sphincterotome, Short Wire	180cm	30mm	Monofilament	5/box

Cannulae




ProForma® HF 4.5 Cannula

Sterile, Disposable

Catalog Number	Tip	Tip O.D.	Catheter O.D.	Length	Guidewire Compatibility	Unit
Precurved						
050907	Curved Taper Tip	4.5 F/1.5mm	7-5 F/1.7mm	190cm	.035"/.89mm	5/box
050906	Straight Taper Tip	4.5 F/1.5mm	7-5 F/1.7mm	190cm	.035"/.89mm	5/box

■ Guidewires




FXWire® Advanced Measurement Guidewire

Sterile, Disposable

Catalog Number	Description	Guidewire O.D. Inches	Tip Length	Wire Length	Tip Type	Torque Vise	Unit
FX35450S	Standard Wire	.035"	5cm	450cm	Straight Tip	No	3/box
FX35450A	Standard Wire	.035"	5cm	450cm	Torquable Angled Tip	Yes	3/box
FX35260S	Standard Wire	.035"	5cm	260cm	Straight Tip	No	3/box
FX35260A	Standard Wire	.035"	5cm	260cm	Torquable Angled Tip	Yes	3/box

■ Guidewires




XWire® Guidewire

Sterile, Disposable

Catalog Number	Description	Guidewire O.D. Inches	Tip Length	Wire Length	Tip Type	Torque Vise	Unit
035450S	Standard Wire	.035"	5cm	450cm	Straight Tip	No	3/box
035450A	Standard Wire	.035"	5cm	450cm	Torquable Angled Tip	Yes	3/box
035260S	Standard Wire	.035"	5cm	260cm	Straight Tip	No	3/box
035260A	Standard Wire	.035"	5cm	260cm	Torquable Angled Tip	Yes	3/box
035450SS	Stiff Wire	.035"	5cm	450cm	Straight Tip	No	3/box
035450AS	Stiff Wire	.035"	5cm	450cm	Torquable Angled Tip	Yes	3/box
035260SS	Stiff Wire	.035"	5cm	260cm	Straight Tip	No	3/box
035260AS	Stiff Wire	.035"	5cm	260cm	Torquable Angled Tip	Yes	3/box
025450S	Standard Wire	.025"	5cm	450cm	Straight Tip	No	3/box
025260S	Standard Wire	.025"	5cm	260cm	Straight Tip	No	3/box

Dilation Balloons



Eliminator® PET Biliary Balloon Dilators

Sterile, Disposable

Catalog Number	Description	Balloon Inflation		Catheter		Guidewire	Unit
		Pressure	Length	O.D.	Length	Compatibility	
000347	4mm x 2cm	150 psi/10.2 atm	2cm	5 F/1.7mm	180cm	.035"/.89mm	1/box
000348	6mm x 2cm	150 psi/10.2 atm	2cm	5 F/1.7mm	180cm	.035"/.89mm	1/box
000349	8mm x 3cm	120 psi/8.2 atm	3cm	5 F/1.7mm	180cm	.035"/.89mm	1/box

Stone Removal Balloons




- One Balloon for Any Size Duct
 - Expandable from 9mm-22mm
- Designed to Reduce Resistance and Bursting
- Dual-Injection Ports
- Proximal and distal injection ports on each device
- Unlike current market-leading biliary stone balloons, FIT® is not made with natural rubber latex.
- CONMED Biliary Assurance Program

FIT® Biliary Stone Balloon

Sterile, Disposable

Not made with natural rubber latex.

Catalog Number	Wire Exchange	Balloon		Catheter		Guidewire	Unit
		Size	Max Volume	O.D.	Length	Compatibility	
Standard							
F922SW	Short	9-22mm	4.5mL	7 F/2.3mm	210cm	.035"/.89mm	1/box
F922LW	Long	9-22mm	4.5mL	7 F/2.3mm	210cm	.035"/.89mm	1/box



DuraGlide™ Triple Lumen Stone Removal Balloons


Long Wire Only

Latex

Sterile, Disposable

	Balloon		Catheter		Guidewire	
Catalog Number	O.D.	Volume	O.D.	Length	Compatability	Unit
Injection Above Balloon						
007123	8.5mm	.75 cc	7-5 F/2.3-1.7mm	200cm	.035"/.89mm	1/box
007124	11.5mm	1.5 cc	7-5 F/2.3-1.7mm	200cm	.035"/.89mm	1/box
007125	15.0mm	3.5 cc	7-5 F/2.3-1.7mm	200cm	.035"/.89mm	1/box
Injection Below Balloon						
007143	8.5mm	.75 cc	7-5 F/2.3-1.7mm	200cm	.035"/.89mm	1/box
007144	11.5mm	1.5 cc	7-5 F/2.3-1.7mm	200cm	.035"/.89mm	1/box
007145	15.0mm	3.5 cc	7-5 F/2.3-1.7mm	200cm	.035"/.89mm	1/box

Fine Needle Aspiration



ClearView™ Endoscopic Ultrasound Needle

Sterile, Disposable

Catalog Number	Description	Gauge	Stylet	Needle	Units
				Diameter	
CV2119R	ClearView™ EUS Needle	19	Round	1.1mm x 2.1mm	5/case
CV1822EB	ClearView™ EUS Needle	22	Extended Bevel	0.7mm x 1.8mm	5/case
CV2722EB	ClearView™ EUS Needle	22	Extended Bevel w/Sheath Stabilizer	0.7mm x 2.7mm	5/case
CV1822R	ClearView™ EUS Needle	22	Round	0.7mm x 1.8mm	5/case
CV1825R	ClearView™ EUS Needle	25	Round	0.5mm x 1.8mm	5/case

All ClearView™ Needles include lockable aspiration syringe and valve

PET Balloon Dilators



Eliminator® PET Balloon Dilators

Disposable

Catalog Number	Description	Inflated OD		Balloon	Catheter	Unit
		mm	French	Length	Length	
000340	Esophageal Balloon Dilator	6mm	18 F	8cm	180cm	1/box
000341	Esophageal Balloon Dilator	8mm	24 F	8cm	180cm	1/box
000342	Esophageal Balloon Dilator	10mm	30 F	8cm	180cm	1/box
000343	Esophageal Balloon Dilator	12mm	36 F	8cm	180cm	1/box
000344	Esophageal Balloon Dilator	15mm	45 F	8cm	180cm	1/box
000345	Esophageal Balloon Dilator	18mm	54 F	8cm	180cm	1/box
000380	Esophageal Balloon Dilator	20mm	60 F	8cm	180cm	1/box
000735	Esophageal Advantage Pack (includes 5 ea. of 15mm)			8cm	180cm	5/box
000736	Esophageal Advantage Pack (includes 5 ea. of 18mm)			8cm	180cm	5/box
000379	Esophageal Advantage Pack (includes 5 ea. of 20mm)			8cm	180cm	5/box
000844	Pyloric/Colonic Balloon	6mm	18 F	4cm	230cm	1/box
000845	Pyloric/Colonic Balloon	8mm	24 F	4cm	230cm	1/box
000846	Pyloric/Colonic Balloon	10mm	30 F	4cm	230cm	1/box
000847	Pyloric/Colonic Balloon	12mm	36 F	4cm	230cm	1/box
000848	Pyloric/Colonic Balloon	15mm	45 F	4cm	230cm	1/box
000849	Pyloric/Colonic Balloon	18mm	54 F	4cm	230cm	1/box
000861	Pyloric/Colonic Balloon	20mm	60 F	4cm	230cm	1/box
000858	Pyloric/Colonic Advantage Pack (includes 5 ea. of 15mm)			4cm	230cm	5/box
000859	Pyloric/Colonic Advantage Pack (includes 5 ea. of 18mm)			4cm	230cm	5/box
000867	Pyloric/Colonic Advantage Pack (includes 5 ea. of 20mm)			4cm	230cm	5/box

PET Balloon Dilators



CONMED Balloon Inflation Device

Catalog Number	Description	Unit
100851	CONMED Balloon Inflation Device	1/box

Multi-Stage Inflation Pressures

For Esophageal and Pyloric/Colonic Balloon Dilators

Balloon Size	Stage One	Stage Two	Stage Three*
6mm/18 F	60 psi/4.1 atm	90 psi/6.1 atm	120 psi/8.2 atm
8mm/24 F	60 psi/4.1 atm	90 psi/6.1 atm	120 psi/8.2 atm
10mm/30 F	45 psi/3.1 atm	68 psi/4.0 atm	90 psi/6.1 atm
12mm/36 F	45 psi/3.1 atm	68 psi/4.0 atm	90 psi/6.1 atm
15mm/45 F	30 psi/2.0 atm	45 psi/3.1 atm	60 psi/4.1 atm
18mm/54 F	25 psi/1.7 atm	38 psi/2.6 atm	50 psi/3.4 atm
20mm/60 F	25 psi/1.7 atm	38 psi/2.6 atm	50 psi/3.4 atm

*Maximum recommended inflation pressure

■ Polyvinyl Dilators

CleanGuide™ Disposable OTW
Esophageal Dilators



CleanGuide™ Disposable OTW Esophageal Dilators

Sterile, Disposable

Catalog Number	Description	Size	Unit
CG-05-3	CleanGuide™ Disposable OTW Esophageal Dilator	5mm	3/cs
CG-06-3	CleanGuide™ Disposable OTW Esophageal Dilator	6mm	3/cs
CG-07-3	CleanGuide™ Disposable OTW Esophageal Dilator	7mm	3/cs
CG-08-3	CleanGuide™ Disposable OTW Esophageal Dilator	8mm	3/cs
CG-09-3	CleanGuide™ Disposable OTW Esophageal Dilator	9mm	3/cs
CG-10-3	CleanGuide™ Disposable OTW Esophageal Dilator	10mm	3/cs
CG-11-3	CleanGuide™ Disposable OTW Esophageal Dilator	11mm	3/cs
CG-12-3	CleanGuide™ Disposable OTW Esophageal Dilator	12mm	3/cs
CG-13-3	CleanGuide™ Disposable OTW Esophageal Dilator	13mm	3/cs
CG-14-3	CleanGuide™ Disposable OTW Esophageal Dilator	14mm	3/cs
CG-15-3	CleanGuide™ Disposable OTW Esophageal Dilator	15mm	3/cs
CG-16-3	CleanGuide™ Disposable OTW Esophageal Dilator	16mm	3/cs
CG-17-3	CleanGuide™ Disposable OTW Esophageal Dilator	17mm	3/cs
CG-18-3	CleanGuide™ Disposable OTW Esophageal Dilator	18mm	3/cs
CG-19-3	CleanGuide™ Disposable OTW Esophageal Dilator	19mm	3/cs
CG-20-3	CleanGuide™ Disposable OTW Esophageal Dilator	20mm	3/cs

■ **CleanGuide™ Disposable Spring Tip Guidewire**

CleanGuide™ Disposable Spring Tip Guidewire




The NEW Sterile CleanGuide™ Disposable Marked Spring Tip Guidewire helps reduce the risk of potential cross contamination that could occur from using a reusable guidewire.



Accessories for Esophageal Dilation

DIS150	CleanGuide™ Disposable Marked Spring Tip Guidewire (Single use)	10/box
--------	---	--------

■ Hemostasis Clips




DuraClip® Hemostasis Clip

Sterile, Disposable

Catalog Number	Description	Sheath Length	Min Channel Diameter	Jaw Opening	Qty/Case
DC0235	DuraClip® for Colonoscope	235cm	2.8mm	11mm	10 per box
DC0235W	DuraClip® for Colonoscope	235cm	2.8mm	16mm	10 per box

■ Hemostasis Probes



BiCap® Superconductor™ Multi-electrode Bipolar Probes

Sterile, Disposable

Catalog Number	Description	Shaft Size	Shaft Length	Tip Length	Color	Unit
BP-5200A	Bipolar, 6 Silver Electrode, Hemostasis Probe	5 Fr	200cm	6.5mm	Blue	1/box
BP-7300A	Bipolar, 6 Silver Electrode, Hemostasis Probe	7 Fr	300cm	6.5mm	Yellow	1/box
BP-7350A	Bipolar, 6 Silver Electrode, Hemostasis Probe	7 Fr	350cm	6.5mm	Yellow	1/box
BP-7350S	Bipolar, 6 Silver Electrode, Hemostasis Probe	7 Fr	350cm	3.3mm	Yellow	1/box
BP-10300A	Bipolar, 6 Silver Electrode, Hemostasis Probe	10 Fr	300cm	6.5mm	Orange	1/box

■ Hemostasis Probes



Beamer AVEO™ SideFire Probe



Beamer AVEO™ Argon Probes

Sterile, Disposable

Catalog Number	Description	Unit
FIA-AVE0-1	Filter Integrated Argon Probe 1.8mm X 160cm	10/box
FIA-AVE0-2	Filter Integrated Argon Probe 1.8mm X 320cm	10/box
FIA-AVE0-3	Filter Integrated Argon Probe 2.3mm X 230cm	10/box
FIA-AVE0-4	Filter Integrated Argon Probe 2.3mm X 320cm	10/box
FIA-AVE0-5	Filter Integrated Argon Probe 3.2mm X 230cm	10/box
FIS-AVE0-3	Filter Integrated Argon SideFire Probe 2.3mm X 230cm	10/box



Beamer AVEO™ Argon Probe

■ Disposable Injection Needles



Click-Tip™ Injection Needle

Sterile, Disposable

Catalog Number	Needle Diameter	Needle Length	Catheter OD	Catheter Length	Unit
04-23-230	19 Gauge Injection Needle	6mm	2.3mm	230cm	10
03-19-230	22 Gauge Injection Needle	6mm	1.9mm	230cm	10
05-19-230	25 Gauge Injection Needle	4mm	1.9mm	230cm	10
06-19-230	25 Gauge Injection Needle	6mm	1.9mm	230cm	10




Flexitip™ Disposable Sclerotherapy Needles

Sterile, Disposable

Catalog Number	Description	Needle		Sheath		Unit
		Length	Gauge	Length	O.D.	
000215	Optic Yellow Tip Needle	4mm	25 g	160cm	2.3mm	5/box
000216	Optic Yellow Tip Needle	5mm	25 g	160cm	2.3mm	5/box
000218	Optic Yellow Tip Colonoscope Needle	5mm	25 g	230cm	2.3mm	5/box

■ Ligation Devices



Stiegmann-Goff™ Endoscopic Esophageal Ligator Kit

Latex

Non-Sterile, Disposable

Catalog Number	Description	Inner Cylinders	Maximum Scope Length	Scope Tip Diameter	Minimum Working Channel O.D.	Unit
200220	Endoscopic Ligator	5	130cm	9 -11mm	2.8mm	5 kits/box
000221	Clearvue™ Endoscopic Ligator	5	130cm	9 -11mm	2.8mm	5 kits/box
000227	Clearvue™ Endoscopic Ligator	10	130cm	9 -11mm	2.8mm	5 kits/box

Each Stiegmann-Goff™ Endoscopic Ligator Includes:

- 1 ea. Monofilament trip wire
- 1 ea. Friction fit adaptors for 9-11mm and 11-13mm scopes
- 5 ea. (200220/000221) inner cylinders with preloaded bands or
- 10 ea. (200221/000227) inner cylinders with preloaded bands

■ Ligation Devices



Stiegmann-Goff™ Bandito™ Hemorrhoid Endoscopic Ligator

Latex

Non-Sterile, Disposable

Catalog Number	Description	Endoscope Compatibility	Maximum Scope Length	Scope Tip Diameter	Minimum Working Channel	Unit
100225	Endoscopic Hemorrhoid Ligator	Colonoscope	210cm	13 -15mm	2.8mm	5 kits/box

- 1 ea. Monofilament trip wire
- 1 ea. Friction fit adaptor for 13-15mm scopes
- 3 ea. Clear inner cylinder with preloaded bands

■ Disposable Biopsy Forceps



Oval Cup



Oval Cup with Needle



Alligator Cup



Alligator Cup with Needle

Precisor® EXL™ Coated Disposable Biopsy Forceps

Featuring Channel-Guard® Low Friction Sheath

Sterile, Coated

Catalog Number	Jaw Configuration	Jaw	Shaft		
		O.D.	Length	Color	Unit
Standard GI					
000381	Alligator Cup	2.3mm	160cm	Yellow	10/box
000382	Alligator Cup with Needle	2.3mm	160cm	Yellow	10/box
000383	Oval Cup	2.3mm	160cm	Yellow	10/box
000384	Oval Cup with Needle	2.3mm	160cm	Yellow	10/box
000385	Alligator Cup	2.3mm	230cm	Blue	10/box
000386	Alligator Cup with Needle	2.3mm	230cm	Blue	10/box
000389	Oval Cup	2.3mm	230cm	Blue	10/box
000388	Oval Cup with Needle	2.3mm	230cm	Blue	10/box
Jumbo GI					
00821C	Oval Cup with Needle	3.1mm	230cm	Blue	10/box
00823C	Alligator Cup with Needle	3.1mm	230cm	Blue	10/box



Oval Cup



Oval Cup with Needle



Alligator Cup



Alligator Cup with Needle

Optibite® Disposable Biopsy Forceps

Sterile, Coated

Catalog Number	Jaw Configuration	Jaw Size	Shaft Length	Color	Unit/Box
100575	Alligator cup w/needle	2.3mm	230cm	Blue	10/box
100576	Alligator cup w/o needle	2.3mm	230cm	Blue	10/box
100577	Oval cup w/needle	2.3mm	230cm	Blue	10/box
100578	Oval cup w/o needle	2.3mm	230cm	Blue	10/box

■ Disposable Biopsy Forceps



Oval Cup



Alligator Cup

Pediatric GI

Sterile, Coated

Catalog Number	Jaw Configuration	Jaw	Shaft			Unit
		O.D.	Length	Color		
100512	Oval Cup	1.8mm	160cm	Blue		10/box
100513	Alligator Cup	1.8mm	160cm	Blue		10/box



Oval Cup



Alligator Cup

Precisor® Hot Disposable Biopsy Forceps

Sterile, Coated

Catalog Number	Configuration	Jaw	Shaft				Unit
		O.D.	Length	Coating	Cord		
000852	Oval Cup	2.3mm	230cm	Red	Olympus®		10/box
000854	Alligator Cup	2.3mm	230cm	Red	Olympus		10/box



Oval Cup



Alligator Cup

Precisor® Broncho® Disposable Biopsy Forceps

Sterile, Disposable

Catalog Number	Description	Shaft Length	Jaw O.D.	Unit
100502	Oval Cup Forceps/Channel-Guard™ Coated Shaft	115cm	1.8mm	10/box
100503	Alligator Cup Forceps/Channel-Guard™ Coated Shaft	115cm	1.8mm	10/box
100506	Oval Cup Forceps	115cm	1.8mm	10/box
100507	Alligator Cup Forceps	115cm	1.8mm	10/box


■ Disposable GI Cytology Brushes



CONMED Cytology Brushes

Sterile, Disposable

Catalog Number	Description	Scope	Brush Diameter	Sheath O.D.	Length	Unit
000104	Sheathed	Multi-Purpose	1.7mm	1.7mm	160cm	20/box
000108	Sheathed	Pediatric Gastroscope	3.0mm	1.7mm	160cm	20/box
000109	Sheathed	Duodenoscope	3.0mm	1.8mm	200cm	20/box
000110	Sheathed	Colonoscope	3.0mm	2.1mm	230cm	20/box
000196	Sheathed	Adult Gastroscope	3.0mm	2.0mm	160cm	20/box




Cytology Brushes

Sterile, Disposable

Catalog Number	Description	Brush	Sheath		Unit
		Diameter	Length	Diameter	
151	Straight Handle	2.5mm	160cm	1.7mm	20/box
151R	Ring Handle	2.5mm	160cm	1.7mm	20/box
163	Straight Handle	3.0mm	160cm	2.4mm	20/box
158	Straight Handle	3.0mm	230cm	2.4mm	20/box
158R	Ring Handle	3.0mm	230cm	2.4mm	20/box

■ Disposable Snares




Beamer® Argon Snare Probe

Sterile, Disposable

Catalog Number	Description	Unit
Beamsnare-1	Beamer® Snare Probe 2.5mm x 160 (15mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-2	Beamer® Snare Probe 2.5mm x 230 (15mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-3	Beamer® Snare Probe 2.5mm x 230 (30mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-4	Beamer® Snare Probe 2.5mm x 230 (15mm, Oval, Monofilament Firm Loop)	5/case
Beamsnare-5	Beamer® Snare Probe 2.5mm x 230 (30mm, Oval, Monofilament Firm Loop)	5/case

■ Disposable Snares




Singular® Polypectomy Snares

Sterile, Disposable

Catalog Number	Loop		Active Cord	Sheath		
	Size	Shape		O.D.	Length	Unit
Firm Loop Wire						
000957	Medium	Crescent	Olympus	2.3mm	230cm	10/box
000959	Medium	Hexagonal	Olympus	2.3mm	230cm	10/box
000963	Large	Oval	Olympus	2.3mm	230cm	10/box
000965	Medium	Oval	Olympus	2.3mm	230cm	10/box
000978	Small (Mini)	Oval	Olympus	2.3mm	230cm	10/box
000991	X-Small	Oval	Olympus	2.3mm	230cm	10/box
Soft Loop Wire						
000977	Large	Oval	Olympus	2.3mm	230cm	10/box
000983	Medium	Oval	Olympus	2.3mm	230cm	10/box
000985	Small (Mini)	Oval	Olympus	2.3mm	230cm	10/box
000987	X-Small	Oval	Olympus	2.3mm	230cm	10/box

■ Disposable Snares




Optimizer® Polypectomy Snares

Sterile, Disposable

Catalog Number	Loop		Active Cord	Sheath		Unit
	Size	Shape		O.D.	Length	
Firm Loop Wire						
000472	Large	Oval	Olympus	2.3mm	230cm	10/box
000474	Medium	Oval	Olympus	2.3mm	230cm	10/box
000476	Small (Mini)	Oval	Olympus	2.3mm	230cm	10/box
000478	X-Small	Oval	Olympus	2.3mm	230cm	10/box
Soft Loop Wire						
000460	Large	Oval	Olympus	2.3mm	230cm	10/box
000462	Medium	Oval	Olympus	2.3mm	230cm	10/box
000464	Small (Mini)	Oval	Olympus	2.3mm	230cm	10/box
000466	X-Small	Oval	Olympus	2.3mm	230cm	10/box

■ Disposable Snares



Orbit-Snare® Rotatable Polypectomy Snares

Sterile, Disposable

Catalog Number	Loop		Active Cord	Sheath		Wire Type	Unit
	Size	Shape		O.D.	Length		
01-15-230	15mm	Rotatable Oval	Universal	2.5mm	230cm	Firm	10/box
01-25-230	25mm	Rotatable Oval	Universal	2.5mm	230cm	Firm	10/box
01-35-230	35mm	Rotatable Oval	Universal	2.5mm	230cm	Firm	10/box

CONMED Polypectomy Snares

Sterile, Disposable

Catalog Number	Loop		Active Cord	Sheath		Unit
	Size	Shape		O.D.	Length	
Firm Loop Wire						
100600	16mm	Oval	Olympus	2.3mm	230cm	10/box
100601	23mm	Oval	Olympus	2.3mm	230cm	10/box
100602	32mm	Oval	Olympus	2.3mm	230cm	10/box
Soft Loop Wire						
100603	16mm	Oval	Olympus	2.3mm	230cm	10/box
100604	23mm	Oval	Olympus	2.3mm	230cm	10/box
100605	32mm	Oval	Olympus	2.3mm	230cm	10/box



Active Cord

Non-Sterile, Reusable

Catalog Number	Active Cord
SN-UAC	Universal Monopolar Active Cord (Fits Olympus style electrode connector)



■ Specimen Devices



Optimizer® Multi-Chamber Polyp Trap

Non-Sterile, Disposable

Catalog Number	Description	Chambers	Unit
CAP132	Multi-Chamber Polyp Trap	4	20/box
CAP132B	Multi-Chamber Polyp Trap	4	120/box

■ Biopsy Valves / Introducer Bite Blocks



ClearGuard™ Biopsy Valves

Sterile, Disposable

Catalog Number	Description	Unit
000522	Biopsy Valve for Olympus and Fujinon	100/case
000523	Biopsy Valve for Pentax	100/case

Biopsy Channel Accessory Introducer

Non-sterile, Disposable

Catalog Number	Description	Unit
000160	Accessory Introducer	6/box



Scope Saver® Bite Blocks

Non-sterile, Disposable

Catalog Number	Description	Lumen Size	Color	Unit
000429 [°]	Large Scope Saver® Disposable Latex Free Bite Block	20 x 27mm	Green	50/box
100429	Small Scope Saver® Disposable Latex Free Bite Block	14 x 20mm	Purple	25/box

[°] Accommodates up to a 60 F American Dilator

CONMED Bite Blocks

Non-sterile, Disposable

Catalog Number	Description	Lumen Size	Color	Unit
001429 [°]	CONMED Large Bite Block Blue Latex Free	20 x 27mm	Blue	100/box

[°] Accommodates up to a 60 F American Dilator

■ Endoscope Valve-Tracking and Pre-Cleaning System


- Helps ensure safe handling of endoscope valves during cleaning, disinfection & storage.
- Provides traceability of valves and other detachable parts, keeping them together with parent endoscope during reprocessing.
- Meets standards of AAMI, American National Standard ST91; 2015 & SGNA Standards of Infection Prevention in Reprocessing of Flexible Gastrointestinal Endoscopes.



CleanCage™ Endoscope Valve-Tracking System

Non-Sterile, Disposable

Catalog Number	Description	Unit
CC8475	Endoscope Valve Cage	100/box
CC8400	Dispensing Unit	1/box



- All-in-one tray, channel cleanser and scope sponge for immediate post-procedure bedside endoscope cleaning.
- Just add water to simplify and standardize your endoscope pre-cleaning processes.
- Meets SGNA Practice Committee Standards of Infection Control.



CleanTray™ Endoscope Pre-Cleaner

Non-Sterile, Disposable

Catalog Number	Description	Unit
CT7357	Bedside Pre-Cleaning Endoscope Kit	100/box
CT7300	Dispensing Unit	4/box

■ CleanGuard™ Disposable Endoscope Valves




CleanGuard™ Disposable Endoscope Valves

Sterile, Disposable

Catalog Number	Description	Unit
200305	Air/Water, Suction and Biopsy Valves for Olympus 140/160/180/190/240/260/290 series Endoscopes	50/case
200310	Air/Water, Suction, Biopsy Valves and Aux Water Jet 140/160/180/190/240/260/290 series Endoscopes	50/case

Cleaning Brushes



Channel Master™ Endoscope Cleaning Brushes

Latex **FREE**

Non-Sterile, Disposable – Biopsy Channel Range 2.0-4.2mm, Shaft O.D. 1.8mm

Catalog Number	Scope	Brush		Shaft Length	Unit
		Diameter	Length		
000610	Pediatric Brush	6.4mm	7.4mm	160cm	50/box
000615	CONMED Double Ended Cleaning Brush	6.4mm	7.4mm	230cm	50/box
000617	CONMED Combination Valve Cleaning Brush	7.9mm/6.4mm	29.2mm/7.4mm	230cm	50/box
000619	CONMED Large Combination Valve Cleaning Brush	9.7mm/6.4mm	45.2mm/7.4mm	230cm	50/box
000639	CONMED Single Ended Cleaning Brush	6.4mm	7.4mm	230cm	50/box



- Designed to clean the elevator region of a duodenoscope
- Unique lavage feature
- Removes up to 47% more soil¹ from the distal tip of your duodenoscope compared to a dual-ended cleaning brush.

Whistle® Endoscope Cleaning Brush

Non-Sterile, Disposable – Biopsy Channel Minimum 2.8mm, Shaft O.D. 1.8mm

Catalog Number	Scope	Brush		Shaft Length	Unit
		Diameter	Length		
000101	Whistle® Cleaning Brush	6.4mm	13.1mm	230cm	50/box

Contains female luer adaptor, hollow shaft, and port hole for fluid injection to soften channel obstructions.

¹ Data on file, PDD1763799



■ Transbronchial Aspiration Needles



WANG® Cytology Needles

Sterile, Disposable

Catalog Number	Biopsy Site	Needle		Unit
		Gauge	Length	
MW-122	Central, Carinal and Mediastinal Regions	22 g	13mm	4/box
MW-222	Central or Peripheral Regions	22 g	13mm	4/box
MW-322	Central, Mediastinal and Hilar Regions	22 g	13mm	4/box
MW-522	Peripheral Region	22 g	13mm	4/box
SW-121	Central, Carinal and Mediastinal Regions	21 g	15mm	4/box
SW-221	Central or Peripheral Regions	21 g	15mm	4/box
SW-521	Peripheral Region	21 g	15mm	4/box
W-220	Central or Peripheral Regions	20 g	15mm	4/box

All WANG® needles have a 1.9mm working diameter and a 130cm working length.



WANG® Histology Needles

Sterile, Disposable

Catalog Number	Biopsy Site	Needle		Unit
		Gauge	Length	
MW-319	Central and Mediastinal Regions	inner - 21 g	3mm	4/box
		outer - 19 g	15mm	
MWF-319	Peripheral Region	19 g	15mm	4/box

All WANG® needles have a 1.9mm working diameter and a 130cm working length.

■ Biopsy Forceps



Oval Cup



Alligator Cup

Precisor Broncho® Disposable Biopsy Forceps

Sterile, Disposable

Catalog Number	Description	Shaft Length	Jaw O.D.	Unit
100502	Oval Cup Forceps/Channel-Guard™ Coated Shaft	115cm	1.8mm	10/box
100503	Alligator Cup Forceps/Channel-Guard™ Coated Shaft	115cm	1.8mm	10/box
100506	Oval Cup Forceps	115cm	1.8mm	10/box
100507	Alligator Cup Forceps	115cm	1.8mm	10/box

Brushes



Cytology Brushes

Sterile, Disposable

Catalog Number	Description	Brush Diameter	Length	Sheath		Unit
				Diameter		
133	Straight Handle	1mm	120cm	1.7mm		20/box
133R	Ring Handle	1mm	120cm	1.7mm		20/box
129	Straight Handle	2mm	120cm	1.7mm		20/box
129R	Ring Handle	2mm	120cm	1.7mm		20/box
149	Straight Handle	3mm	120cm	1.7mm		20/box
149R	Ring Handle	3mm	120cm	1.7mm		20/box



Catalog Number	Description	Scope	Brush Diameter	Length	Sheath		Unit
					Diameter		
000104	Sheathed	Multi-Purpose	1.7mm	160cm	1.7mm		20/box
000107	Sheathed	Bronchoscope	3.0mm	115cm	1.7mm		20/box
000139	Harrell™ (unsheathed)	Bronchoscope	4.0mm	110cm	1.7mm		20/box
000140	Harrell™ (unsheathed)	Bronchoscope	7.0mm	110cm	1.2mm		20/box

Brushes



NB-120




PB-120

Needle Tip and Peripheral Cytology Brushes

Sterile, Disposable

Catalog Number	Description	Brush Diameter	Length	Sheath		Unit
				Diameter	Length	
NB-120*	Needle Tip Cytology Brush	1.7mm	11mm	1.9mm	140cm	4/box
PB-120°	Peripheral Radiopaque Tip Cytology Brush	1.7mm	11mm	1.9mm	140cm	4/box

* Needle tip cytology brush is used for retrieving cytology specimens from endobronchial lesions, peripheral lung nodules, or lung masses.
° The metal tipped sheath of the peripheral brush allows for flourosopic visualization during cytologic collection.



Microbiology Brush

Sterile, Disposable

Catalog Number	Description	Brush Diameter	Length	Sheath		Unit
				Diameter	Length	
130	Microbiology Brush	1mm	110cm	1.8mm		5/box


Designed to obtain contamination free specimens, the brush is protected behind a wax plug which is expelled before sampling.

Cleaning Brushes

Channel Master™ Endoscope Cleaning Brush

Non-Sterile, Disposable

Catalog Number	Scope	Biopsy Channel Range	Shaft Length	Unit
000610	Bronchoscope	2.0 - 4.2mm	160cm	50/box *



Reusable Cleaning Brush

Non-Sterile, Reusable

Catalog Number	Scope	Brush		Shaft Length	Unit
		Diameter	Length		
380	Bronchoscope	3mm	9.5mm	115cm	1/box

The reusable scope cleaning brush features a controllable ring handle and a Teflon® coated spring sheath.

Generators & Argon Probes



Beamer AVEO™

Catalog Number	Description
AVEO-ESU	Electrosurgical Generator
AVEO-ABC	Argon Module
AVEO-1W	Wired Foot Pedal
AVEO-CART-PKG	Cart (no tank holder)
AVEO-TANK-LRG	Large Tank Holder (8.25"/20.9cm)
CB200-A02	Argon Pressure Reducer




Beamer AVEO™ Argon Probes

Sterile, Disposable

Catalog Number	Description	Unit
FIA-AVEO-1	Filter Integrated Argon Probe 1.8mm X 160cm	10/box
FIA-AVEO-2	Filter Integrated Argon Probe 1.8mm X 320cm	10/box
FIA-AVEO-3	Filter Integrated Argon Probe 2.3mm X 230cm	10/box
FIA-AVEO-4	Filter Integrated Argon Probe 2.3mm X 320cm	10/box
FIA-AVEO-5	Filter Integrated Argon Probe 3.2mm X 230cm	10/box
FIS-AVEO-3	Filter Integrated Argon SideFire Probe 2.3mm X 230cm	10/box

Beamer® System CE600

Catalog Number	Description
CE200	ESU
CB200	Argon Module
CC200	Cart
CE200-A05	Beamer® Double Foot Pedal with Swap-Mode Button and Foot Lift Bar
CB200-A02	Argon Pressure Reducer
CB200-A01	Beamer® Connecting Cable for Argon Probes
CE200-A04	Connecting Cable for Reusable Return Electrodes
SN-UAC	Universal Monopolar Active Cord



Generators & Argon Probes



Beamer® Argon Probes

Sterile, Disposable


Catalog Number	Description	Unit
A-Beam-1	Beamer® Probe 1.8mm x 160cm	10/case
A-Beam-2	Beamer® Probe 1.8mm x 320cm	10/case
A-Beam-3	Beamer® Probe 2.3mm x 230cm	10/case
A-Beam-4	Beamer® Probe 2.3mm x 320cm	10/case
A-Beam-5	Beamer® Probe 3.2mm x 230cm	10/case



Beamer® Argon SideFire Probe

Sterile, Disposable

Catalog Number	Description	Unit
S-Beam-3	Beamer® Argon SideFire Probe 2.3mm x 230cm	10/case



Beamer® Argon Snare Probes

Sterile, Disposable

Catalog Number	Description	Unit
Beamsnare-1	Beamer® Snare Probe 2.5mm x 160 (15mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-2	Beamer® Snare Probe 2.5mm x 230 (15mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-3	Beamer® Snare Probe 2.5mm x 230 (30mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-4	Beamer® Snare Probe 2.5mm x 230 (15mm, Oval, Monofilament Firm Loop)	5/case
Beamsnare-5	Beamer® Snare Probe 2.5mm x 230 (30mm, Oval, Monofilament Firm Loop)	5/case

Generators & Argon Probes



BiCap® III Electrosurgical Generator

Catalog Number	Description
60-8512-SYS	BiCap® III System with Monopolar and Bipolar Footswitches*
60-8512-120	BiCap® III Electrosurgical Generator
60-6700-001	Monopolar Footswitch
60-5103-001	BiPolar Footswitch
60-8510-001	BiCap® III Pump Mounting Bracket

*Includes 60-8512-120, 60-6700-001 and 60-5103-001.


Beamer AVEO™ Irrigation Pump*

Catalog Number	Description
AVEOPUMP	Beamer AVEO™ Irrigation Pump
PUMP-1	Beamer AVEO™ Irrigation Pump Foot Pedal Replacement Part
PUMP-2	Beamer AVEO™ Irrigation Pump Head Replacement Part
PUMP-3	Beamer AVEO™ Irrigation Pump Bottle Holder (US/EU) Replacement Part

*Not available in EU and Australia



Generators & Argon Probes



Coaxial Bipolar Adaptors

Non-Sterile, Reusable

Catalog Number	Description	Unit
133001	Coaxial Bipolar Adapter (For use with the System 7500® & System 7550®)	1/box
133002	Coaxial Bipolar Adapter * (For use with the System 5000™ & System 2450®)	1/box

* For use with Everest™ and ValleyLab™ Generators



Active Cord

Non-Sterile, Reusable

Catalog Number	Description	Unit
SN-UAC	Universal Monopolar Active Cord (Fits Olympus style electrode connector)	1/box



ABCFlex™ Argon Probes

Sterile, Disposable

Catalog Number	Description	Unit
133023	ABCFlex™ Probe 2.3mm (7Fr.) x 220cm, Original	10/case
133270	ABCFlex™ Probe 2.3mm (7Fr.) x 270cm, Original	10/case
133327	ABCFlex™ Probe 2.3mm (7Fr.) x 270cm, System 7550®	10/case

Please be advised that the following devices or components provided by CONMED contain natural rubber latex, which may cause allergic reactions in some individuals:

Product	Catalog Number	Device/Component
Stiegmann-Goff™ Endoscopic Ligator Kit	000221, 000227, 200220-221	"O" ring contains latex
Stiegmann-Goff™ Bandito™ Endoscopic Ligator	100225	"O" ring contains latex
DuraGlide3™ Stone Removal Balloon	007123-007125, 007143-007145	Latex balloon
DuraGlide™ Stone Removal Balloon	000635, 000637-638, 000646-647	Latex balloon
Disposable Percutaneous Biopsy Instruments	000402, 000439, 000441	Syringe Plunger

<u>Catalog Number</u>	<u>Page</u>	<u>Catalog Number</u>	<u>Page</u>	<u>Catalog Number</u>	<u>Page</u>
000104	22	000474	25	007125	6
000107	36	000476	25	007143	6
000108	22	000478	25	007144	6
000109	22	000522	29	007145	6
000110	22	000523	29	00821C	18
000139	36	000610	32	00823C	18
000140	36	000615	32	01-15-230	26
000160	29	000617	32	01-25-230	26
000196	22	000619	32	01-35-230	26
000215	17	000639	32	025260S	4
000216	17	000735	10	025450S	4
000218	17	000736	10	03-19-230	16
000221	18	000844	10	035260A	5
000227	18	000845	10	035260AS	5
000340	10	000846	10	035260S	5
000341	10	000847	10	035260SS	6
000342	11	000848	10	035450A	6
000343	11	000849	10	035450AS	6
000344	11	000852	22	035450S	6
000345	10	000854	22	035450SS	6
000347	8	000858	10	04-23-230	16
000348	8	000859	10	05-19-230	16
000349	8	000861	10	050906	3
000379	10	000867	10	050907	3
000380	10	000957	24	06-19-230	18
000381	19	000959	24	100225	18
000382	19	000963	24	100429	29
000383	19	000965	24	100502	23
000384	19	000977	24	100503	23
000385	19	000978	24	100506	23
000386	19	000983	24	100507	23
000388	19	000985	24	100512	23
000389	19	000987	24	100513	23
000429	22	000991	24	100575	18
000460	25	001429	29	100576	18
000462	25	007104	3	100577	18
000464	25	007105	3	100578	18
000466	25	007123	6	100600	26
000472	25	007124	6	100601	26

CATALOG NUMBER INDEX

Catalog Number	Page	Catalog Number	Page	Catalog Number	Page
100602	26	A-Beam-4	40	FX35260A	5
100603	26	A-Beam-5	40	FX35260S	5
100604	26	Beamsnare-1	40	FX35450A	5
100605	26	Beamsnare-2	40	FX35450S	5
100851	10	Beamsnare-3	40		
129	40	Beamsnare-4	40		
129R	40	Beamsnare-5	40		
130	41	BP-10300A	16		
133	40	BP-5200A	16		
133001	18	BP-7300A	16		
133002	18	BP-7350A	16		
133023	18	BP-7350S	16		
133270	18	CAP132	32		
133327	18	CAP132B	32		
133R	40	CB200	39		
149	40	CB200-A01	39		
149R	40	CB200-A02	39		
151	25	CC200	39		
151R	25	CC8475	34		
158	25	CE200	39		
158R	25	CE200-A04	39		
163	25	CE200-A05	39		
200220	20	CT7357	34		
200305	35	CV1822EB	11		
200310	35	CV1822R	11		
60-5103-001	41	CV1825R	11		
60-6700-001	41	CV2119R	11		
60-8512-120	41	CV2722EB	11		
60-8512-SYS	41	DC0235	14		
7104AC	3	DC0235W	14		
7105AC	3	DIS150	15		
AVEO-ABC	41	F922SW	8		
AVEO-CART-PKG	41	F922LW	8		
AVEO-ESU	41	FIA-AVEO-1	16		
AVEO-TANK-LRG	41	FIA-AVEO-2	16		
AVEOPUMP	41	FIA-AVEO-3	16		
A-Beam-1	40	FIA-AVEO-4	16		
A-Beam-2	40	FIA-AVEO-5	16		
A-Beam-3	40	FIS-AVEO-3	16		

CATALOG NUMBER INDEX

Catalog Number	Page
MW-122	38
MW-222	38
MW-319	38
MW-322	38
MW-522	38
MWF-319	38
NB-120	41
PB-120	41
PUMP-1	41
PUMP-2	41
PUMP-3	41
S-BEAM-3	40
SLW20-5096-5.....	4
SLW30-5096-5.....	4
SN-UAC	18
SSW20-50912-5.....	4
SSW30-50912-5.....	4
SW-121	38
SW-221	38
SW-521	38
T2043F	2
T3043F	2
W-220	38



Imagine what we could do together.

CONMED Corporation | 11311 Concept Blvd., Largo, FL 33773 | www.CONMED.com



Service & Support